



WE CARE LIKE FAMILY™

EVER-SOFT™ FOAM WEDGE & COCCYX WHEELCHAIR CUSHIONS



- High-density foam provides comfort and support
- Cushion encourages healing via active evaporation of moisture and cooling at patients' pressure points
- Includes removable, low-shear cover with non-skid bottom
- Convenient, adjustable strap secures to wheelchair
- Weight capacity 300 lb.
- Wedge (10645-10646) helps prevent forward sliding
- Coccyx (10655-10656) has a contoured design to relieve pressure off the coccyx
- 1-Year Warranty
- Meets CAL 117-2013

Reorder No.	Description	Size	Qty./Bx	Bx/Cs
10645	Ever-Soft™ Wedge Foam Cushion	18"×16"×4"×2"	1	
10645-4	Ever-Soft™ Wedge Foam Cushion	18"×16"×4"×2"	1	4
10646	Ever-Soft™ Wedge Foam Cushion	20"×16"×4"×2"	1	
10646-4	Ever-Soft™ Wedge Foam Cushion	20"×16"×4"×2"	1	4
10655	Ever-Soft™ Coccyx Foam Cushion	18"×16"×3"	1	
10655-4	Ever-Soft™ Coccyx Foam Cushion	18"×16"×3"	1	4
10656	Ever-Soft™ Coccyx Foam Cushion	20"×16"×3"	1	
10656-4	Ever-Soft™ Coccyx Foam Cushion	20"×16"×3"	1	4

www.dynarex.com • 888-396-2739

Dynarex Corporation • 10 Glenshaw Street • Orangeburg, NY 10962

R210316



© 2021 Dynarex Corporation