## ELECTRIC PATIENT LIFT Sit-To-Stand Patient Lift

## FEATURES AND BENEFITS

- Durable steel frame construction
- Equipped with hand control, diagnostic display, removable battery pack, and emergency stop button
- Dual 3" front casters and lockable 5" rear casters
- Sling hook-up height adjusts between 33.6"-70"
- Base legs adjust easily into open position with padded shifter handle
- Optional accessories include wall charger and slings
- Soft, adjustable knee pad for lower leg support
- Foot plate with anti-slip surface
- Conveniently located handles for patient to hold during transport





Reorder No.	Weight Capacity	Height at Sling Hook-up		Base Width		Base Height (Clearance)	Base Length	HCPCS
		Max.	Min.	Open	Closed	(Olearance)		
10702	450 lb.	70"	33.6"	40"	25"	4.5"	48"	E0635

## ALSO AVAILABLE:

Sit-to-Stand Slings (10746-10747) and Transfer Slings (10748)

Product Description	SIZ	Quantity	
	Large	X-Large	quantity
Deluxe Sit-to-Stand Sling	10746LG	10746XL	1/Bx
Sit-to-Stand Sling	10747LG	10747XL	1/Bx
Transfer Sling- Solid	10748LG	10748XL	1/Bx



## **dynarex**<sup>®</sup> is your one-stop medical products brand since 1967

For more information contact your local Dynarex sales representative OR Call us at: (888) 396-2739 • www.dynarex.com