



WE CARE LIKE FAMILY™



## Sterile Alcohol Prep Pads

### For Optimal Antibacterial Action

- For topical cleansing prior to injection or venipuncture and manufactured according to current FDA regulations, GMP and USP requirements
- Provides rapid-acting skin antiseptics in a sterile formula
- Sterile alcohol prep pads are saturated with 70% isopropyl alcohol
- Soft and absorbent, non-woven pad
- Gamma sterilized and not made with natural rubber latex



Reorder No.	Size	Length Folded	Length Unfolded	Width	Quantity
1113	Medium	1.18"	2.36"	1.1"	10/200/Cs
1114	Medium	1.18"	2.69"	1.1"	20/100/Cs
1116	Large	1.77"	3.54"	1.18"	10/100/Cs

\*Sizes will vary slightly based on product specification tolerances

HCPCS: A4245



Find us on:



**dynarex**® is your one-stop medical products brand since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • [www.dynarex.com](http://www.dynarex.com)

© 2018 Dynarex Corporation

R180802