

# LONGTERM CARE BEDS

## D100 Series Beds



WE CARE LIKE FAMILY™

### FEATURES AND BENEFITS

- Constructed from welded steel
- Grid deck and mattress retainers accommodate both foam or alternating pressure mattresses
- Ultra-quiet DC motor has 4 functions: head, foot, auto contour, and bed height adjustment
- Deck to floor range adjusts from 10.5" – 24.5"
- Equipped with 4 locking casters and wall bumper guard
- 450 lb. maximum weight capacity

*Warranty: 15 years on all welds, 3 years on electrical and mechanical, 1 year on all other parts and components from original date of purchase.*



### ALSO AVAILABLE:

- **Tool-free, Removable Headboards/Footboards:**  
*Available in Mahogany, Amber Cherry, and Light Oak or Composite Plastic*
- Metal Swing Rail, Countoured Swing Rail, or Assist Bar
- Footboard with Embedded Nurses Control
- Battery Backup



Mahogany



Amber Cherry



Light Oak

Nurses Control Panel



Find Us On:



**dynarex**® is your one-stop medical products brand since 1967

For more information contact your local Dynarex sales representative OR  
Call us at: (888) 396-2739 • [www.dynarex.com](http://www.dynarex.com)