





WE CARE LIKE FAMILY™

Advanced Wound Care Dressing DynaGinate™ AG Silver Calcium Alginate Dressing



Used on infected wounds and wounds suspected of possible infection

- Provides 3 days of continuous antibacterial protection
- Used on infections; pressure ulcers; leg ulcers; diabetic wounds; donor sites; moderate to heavy draining, partial and full thickness wounds
- Same as Calcium Alginate dressing, with antibacterial Silver
- Absorbs 12 times its own weight in fluid

Reorder No.	DynaGinate™ AG Silver Calcium Alginate Dressing Description	Quantity
3030	2" x 2" – Silver Calcium Alginate Dressing	10/10/Cs
3031	4" x 5" – Silver Calcium Alginate Dressing	10/10/Cs
3032	4" x 8" – Silver Calcium Alginate Dressing	10/5/Cs
3033	2 g x 30 cm – Silver Calcium Alginate Rope Dressing	10/5/Cs

3030 HCPCS: A6196, 3031 & 3032 HCPCS: A6197, 3033 HCPCS: A6199

Information provided about Dynarex advanced wound care products should be considered as general guidelines to product usage and none of the information is intended to be construed as medical advice. Medical conditions should be evaluated and treated by licensed practitioners. Please see your local wound care clinician for treatment of your wound.

PLUS+

Many Additional
Advanced Wound
Care Dressings
Available From
Dynarex!



SEE OTHER REFERENCE
MATERIALS ON OUR WEBSITE
FOR MORE INFORMATION



Find us on:



dynarex® is your one-stop medical supply company since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarex.com