



CLOSED SUCTION CATHETERS ENDOTRACHEAL DOUBLE SWIVEL ELBOW • LENGTH: 54.0 CM/21.3 IN



Dynarex Resp-O₂™ Closed Suction Endotracheal Catheters are disposable respiratory devices designed to safely remove secretions from the respiratory tract of ventilator-dependent patients. The device allows ventilation to continue during suctioning, helps prevent cross-contamination, and protects caregivers from airborne pathogens.

- Protects patient from contact with outside pathogens and prevents contamination of the surrounding environment
- Smooth, flexible tubing minimizes patient discomfort
- Pressure wiper seal helps maintain PEEP
- Flex Tube provides adjustment and flexibility to the ventilator circuit
- Designed for 24-hr, single-patient use
- Does not contain DEHP

10 Fr

Double Swivel Elbow

REORDER NO. 31030

ADULT/PEDIATRIC

ACCESSORIES

15 mm x 22 mm Flex Tube
Change of Day

PACKAGE QUANTITY
20/Case

12 Fr

Double Swivel Elbow

REORDER NO. 31031

ADULT/PEDIATRIC

ACCESSORIES

15 mm x 22 mm Flex Tube
Change of Day

PACKAGE QUANTITY
20/Case

14 Fr

Double Swivel Elbow

REORDER NO. 31032

ADULT

ACCESSORIES

15 mm x 22 mm Flex Tube
Change of Day

PACKAGE QUANTITY
20/Case

16 Fr

Double Swivel Elbow

REORDER NO. 31033

ADULT

ACCESSORIES

15 mm x 22 mm Flex Tube
Change of Day

PACKAGE QUANTITY
20/Case



WE CARE LIKE FAMILY™

1.888 DYNAREX

DYNAREX.COM

10 GLENSHAW STREET | ORANGEBURG, NY 10962

Contact your local sales representative or call to find your nearest distributor.



©2021 Dynarex Corporation

R210303