





## **NASAL OXYGEN CANNULAS**

## Cushion Tip with Ultra Soft Ear Foam

Dynarex Resp-O<sub>2</sub>™ Nasal Oxygen Cannulas with Cushion Tips and Ultra Soft Ear Foam are designed to efficiently deliver supplemental oxygen to respiratory patients. The cannulas come with either Standard or Universal connectors and feature light green oxygen supply tubing, designed to help improve tube visibility.

- Adjustable foam ear cushions and lightweight cushioned nasal prongs increase comfort and help minimize the risk of irritation
- Crush-resistant, star lumen oxygen supply tubing
- Individually wrapped
- Not made with natural rubber latex or DEHP

7' (2.1 m) Crush **Resistant Tubing** 

STANDARD CONNECTOR

REORDER NO. 33150

**PEDIATRIC** 

PACKAGE QUANTITY 50/Case

7' (2.1 m) Crush **Resistant Tubing** 

UNIVERSAL CONNECTOR

REORDER NO. 33151

PEDIATRIC

PACKAGE QUANTITY 50/Case

4' (1.2 m) Crush **Resistant Tubing** 

STANDARD CONNECTOR

REORDER NO. 33152

ADULT

PACKAGE QUANTITY 50/Case

7' (2.1 m) Crush **Resistant Tubing** 

STANDARD CONNECTOR

REORDER NO. 33153

**ADULT** 

**PACKAGE QUANTITY** 50/Case

7' (2.1 m) Crush **Resistant Tubing** 

UNIVERSAL CONNECTOR

REORDER NO. 33154

**ADULT** 

PACKAGE QUANTITY 50/Case

14' (4.3 m) Crush **Resistant Tubing** 

STANDARD CONNECTOR

REORDER NO. 33155

**ADULT** 

PACKAGE QUANTITY 50/Case

25' (7.6 m) Crush **Resistant Tubing** 

STANDARD CONNECTOR


**REORDER NO. 33156** 

ADULT

25/Case

PACKAGE QUANTITY











Contact your local sales representative or call to find your nearest distributor





