





Dynarex Resp-O<sub>2</sub>™ Nasal Oxygen Cannulas with Cushion Tips and Advantage Ear Foam are designed to efficiently deliver supplemental oxygen to respiratory patients. These translucent cannulas come with either Standard or Universal connectors

and feature light green oxygen supply tubing, designed to help

Adjustable sponge-foam ear wraps and cushioned nasal prongs provide added comfort and help minimize the risk of

- Crush-resistant, star lumen oxygen supply tubing
- Individually wrapped
- Not made with natural rubber latex or DEHP

7' (2.1 m) Crush **Resistant Tubing** 

STANDARD CONNECTOR

**REORDER NO. 33170** 

**PEDIATRIC** 

PACKAGE QUANTITY 50/Case

7' (2.1 m) Crush **Resistant Tubing** 

UNIVERSAL CONNECTOR

**REORDER NO. 33171** 

PEDIATRIC

PACKAGE QUANTITY 50/Case

ADULT

PACKAGE QUANTITY 50/Case

4' (1.2 m) Crush

**Resistant Tubing** 

STANDARD CONNECTOR

**REORDER NO. 33172** 

7' (2.1 m) Crush **Resistant Tubing** 

STANDARD CONNECTOR

REORDER NO. 33173

**ADULT** 

improve tube visibility.

PACKAGE QUANTITY 50/Case

7' (2.1 m) Crush **Resistant Tubing** 

UNIVERSAL CONNECTOR

**REORDER NO. 33174** 

**ADULT** 

PACKAGE QUANTITY 50/Case

14' (4.3 m) Crush **Resistant Tubing** 

STANDARD CONNECTOR

REORDER NO. 33175

**ADULT** 

PACKAGE QUANTITY 50/Case

25' (7.6 m) Crush **Resistant Tubing** 


STANDARD CONNECTOR

**REORDER NO. 33176** 

**ADULT** 

PACKAGE QUANTITY 25/Case







DYNAREX.COM



Contact your local sales representative or call to find your nearest distributor





