

Equal Length EtCO₂ Sampling Line & O₂ Tubing

EtCO₂ NASAL SAMPLING CANNULAS

Dynarex Resp-O₂™ EtCO₂ Nasal Sampling Cannulas with Cushion Nasal Tips are designed to help monitor exhaled carbon dioxide during respiratory treatment. Each soft, cushioned cannula features a permanent barrier between the nasal prongs so that EtCO₂ can be sampled from one nostril while oxygen continues to flow to the other nostril.

- Permanent barrier between nasal prongs allows for EtCO₂ sampling with uninterrupted oxygen delivery
- Constructed from soft materials for added patient comfort
- Crush-resistant tubing helps prevent occlusion and kinking
- Features EtCO₂ sampling line and O₂ tubing of equal length
- Transparent design allows for clear visualization
- Not made with natural rubber latex or DEHP



Female Luer Connector Male Luer Connector

Standard Tubing Connector

Universal Tubing Connector 7' (2.1 m) Tubing

MALE LUER CONNECTOR STANDARD O2 TUBING CONNECTOR

REORDER NO. 33180

INFANT

PACKAGE QUANTITY 25/Case

7' (2.1 m) Tubing

MALE LUER CONNECTOR STANDARD O2 TUBING CONNECTOR

REORDER NO. 33181

PEDIATRIC

PACKAGE QUANTITY 25/Case

7' (2.1 m) Tubing

MALE LUER CONNECTOR STANDARD O2 TUBING CONNECTOR

REORDER NO. 33182

ADULT

PACKAGE QUANTITY 25/Case

7' (2.1 m) Tubing

FEMALE LUER CONNECTOR STANDARD O2 TUBING CONNECTOR

REORDER NO. 33183

INFANT

PACKAGE QUANTITY

25/Case

7' (2.1 m) Tubing

FEMALE LUER CONNECTOR STANDARD O2 TUBING CONNECTOR

REORDER NO. 33184

PEDIATRIC

PACKAGE QUANTITY 25/Case

7' (2.1 m) Tubing

FEMALE LUER CONNECTOR STANDARD O2 TUBING CONNECTOR

REORDER NO. 33185

ADULT

PACKAGE QUANTITY 25/Case

10' (3.0 m) Tubing

MALE LUER CONNECTOR STANDARD O2 TUBING CONNECTOR

REORDER NO. 33186

PEDIATRIC

PACKAGE QUANTITY 25/Case

10' (3.0 m) Tubing

MALE LUER CONNECTOR STANDARD O2 TUBING CONNECTOR

REORDER NO. 33187

ADULT

PACKAGE QUANTITY 25/Case

14' (4.2 m) Tubing

MALE LUER CONNECTOR STANDARD O2 TUBING CONNECTOR

REORDER NO. 33188

PEDIATRIC

PACKAGE QUANTITY 25/Case

14' (4.2 m) Tubing

MALE LUER CONNECTOR STANDARD O2 TUBING CONNECTOR

REORDER NO. 33189

ADULT

PACKAGE QUANTITY 25/Case

10' (3.0 m) Tubing

FEMALE LUER CONNECTOR STANDARD O2 TUBING CONNECTOR

REORDER NO. 33190

PEDIATRIC

PACKAGE QUANTITY 25/Case

10' (3.0 m) Tubing

FEMALE LUER CONNECTOR STANDARD O2 TUBING CONNECTOR

REORDER NO. 33191

ADULT

PACKAGE QUANTITY 25/Case

14' (4.2 m) Tubing

FEMALE LUER CONNECTOR STANDARD O2 TUBING CONNECTOR

REORDER NO. 33192

PEDIATRIC

PACKAGE QUANTITY 25/Case

14' (4.2 m) Tubing

FEMALE LUER CONNECTOR STANDARD O2 TUBING CONNECTOR

REORDER NO. 33193

ADULT

PACKAGE QUANTITY 25/Case

7' (2.1 m) Tubing

MALE LUER CONNECTOR UNIVERSAL O2 TUBING CONNECTOR

REORDER NO. 33194

PEDIATRIC

PACKAGE QUANTITY 25/Case

7' (2.1 m) Tubing

MALE LUER CONNECTOR UNIVERSAL O2 TUBING CONNECTOR

REORDER NO. 33195

ADULT

PACKAGE QUANTITY 25/Case

7' (2.1 m) Tubing

FEMALE LUER CONNECTOR UNIVERSAL O2 TUBING CONNECTOR

REORDER NO. 33196

ADULT

PACKAGE QUANTITY 25/Case

10' (3.0 m) Tubing

FEMALE LUER CONNECTOR UNIVERSAL O2 TUBING CONNECTOR

REORDER NO. 33197

ADULT

PACKAGE QUANTITY 25/Case

WE CARE LIKE FAMILY

1.888 DYNAREX



DYNAREX.COM



O 10 GLENSHAW STREET | ORANGEBURG, NY 10962

Contact your local







R210426

sales representative or call to find your nearest distributor