




COMFORT-FOAM TRACH TUBE HOLDERS

Dynarex Resp-O₂[™] Comfort-Foam Trach Tube Holders are designed to stabilize and hold tracheostomy tubes securely in place to ensure uninterrupted airflow to trach patients. The Trach Tube Holder fastens securely with hook-and-loop closures and features soft, Comfort-Foam padding to help keep skin dry and promote optimal patient comfort.

- Helps to stabilize trach tube and reduce risk of dislodgement
- Comfort-Foam cushioning helps to minimize skin irritation
- Features high-quality, USA-made materials
- Hook-and-loop closures are easy to adjust
- Not made with natural rubber latex or DEHP



<p>1-Piece</p> <p>SMALL</p> <p>REORDER NO. 35020</p> <p>NEONATAL & INFANT</p> <p>NECK 7"-9"/18-23 cm</p> <p>PACKAGE QUANTITY 5 Boxes/Case • 12/box</p>	<p>1-Piece</p> <p>MEDIUM</p> <p>REORDER NO. 35021</p> <p>ADOLESCENT & ADULT</p> <p>NECK 9"-14"/23-35 cm</p> <p>PACKAGE QUANTITY 5 Boxes/Case • 12/box</p>	<p>1-Piece</p> <p>LARGE</p> <p>REORDER NO. 35022</p> <p>ADULT</p> <p>NECK 13"-19"/33-48 cm</p> <p>PACKAGE QUANTITY 5 Boxes/Case • 12/box</p>	<p>1-Piece</p> <p>X-LARGE</p> <p>REORDER NO. 35023</p> <p>ADULT & BARIATRIC</p> <p>NECK 16-24"/41-61 cm</p> <p>PACKAGE QUANTITY 5 Boxes/Case • 12/box</p>
--	---	--	---

<p>Anti-Disconnect • 1-Piece</p> <p>SMALL</p> <p>REORDER NO. 35024</p> <p>NEONATAL & INFANT</p> <p>NECK 7"-9"/18-23 cm</p> <p>PACKAGE QUANTITY 5 Boxes/Case • 12/box</p>	<p>Anti-Disconnect • 1-Piece</p> <p>MEDIUM</p> <p>REORDER NO. 35025</p> <p>ADOLESCENT & ADULT</p> <p>NECK 9"-14"/23-35 cm</p> <p>PACKAGE QUANTITY 5 Boxes/Case • 12/box</p>	<p>Anti-Disconnect • 1-Piece</p> <p>LARGE</p> <p>REORDER NO. 35026</p> <p>ADULT</p> <p>NECK 13"-19"/33-48 cm</p> <p>PACKAGE QUANTITY 5 Boxes/Case • 12/box</p>	<p>2-Piece</p> <p>SMALL</p> <p>REORDER NO. 35027</p> <p>NEONATAL & INFANT</p> <p>MAXIMUM NECK SIZE 9"/23 cm</p> <p>PACKAGE QUANTITY 5 Boxes/Case • 10/box</p>	<p>2-Piece</p> <p>MEDIUM</p> <p>REORDER NO. 35028</p> <p>PEDIATRIC</p> <p>MAXIMUM NECK SIZE 18"/46 cm</p> <p>PACKAGE QUANTITY 5 Boxes/Case • 10/box</p>	<p>2-Piece</p> <p>LARGE</p> <p>REORDER NO. 35029</p> <p>ADULT</p> <p>MAXIMUM NECK SIZE 19.5"/50 cm</p> <p>PACKAGE QUANTITY 5 Boxes/Case • 10/box</p>	<p>2-Piece</p> <p>X-LARGE</p> <p>REORDER NO. 35030</p> <p>BARIATRIC</p> <p>MAXIMUM NECK SIZE 28"/71 cm</p> <p>PACKAGE QUANTITY 10/Case</p>
--	---	--	---	---	--	--



WE CARE LIKE FAMILY[™]

1.888 DYNAREX

DYNAREX.COM

10 GLENSHAW STREET | ORANGEBURG, NY 10962

Contact your local sales representative or call to find your nearest distributor.



©2021 Dynarex Corporation

R210823