

OXYGEN REGULATORS CGA 870-STYLE

Dynarex Resp-O₂™ Oxygen Flow Regulators are designed to safely and efficiently administer and regulate the flow of oxygen to respiratory patients at adjustable flow-rate levels. They feature an anodized aluminum body with a brass core and brass-encased component parts for optimal safety and reliability. Some select models feature all-brass construction. Designed to fit any oxygen tank with a CGA 870-style fitting.

- Single-selector knob allows for effortless adjustment
- Features easy-to-read flow-rate dial
- Brass core construction provides increased safety and reliability
- Complies to CGA 870 standard
- Not made with natural rubber latex or DEHP



BARB
Outlet

36350-36353



BARB
Outlet

36354-36357



DISS
Outlet

36358-36359



BARB
Outlet

36360



DISS
Outlet

BARB
Outlet

36361

ALSO AVAILABLE:



Cylinder Wrench

REORDER NO. 36600



WE CARE LIKE FAMILY™

1.888 DYNAREX

DYNAREX.COM

10 GLENSHAW STREET | ORANGEBURG, NY 10962

©2021 Dynarex Corporation

Contact your local sales representative or call to find your nearest distributor.



R210823



OXYGEN REGULATORS CGA 870-STYLE FLOW RANGE CHART

Reorder No.	Product Features	Flow Range	Dial settings												Qty.
			1	2	3	4	5	6	7	8	9	10	11	12	
36350	BARB OUTLET • BRASS-SLEEVED	0-4 LPM	0	1/32	1/16	1/8	1/4	1/2	3/4	1	1 1/2	2	3	4	5/CASE
36351	BARB OUTLET • BRASS-SLEEVED	0-8 LPM	0	1/8	1/2	1	1 1/2	2	2 1/2	3	4	5	6	8	5/CASE
36352	BARB OUTLET • BRASS-SLEEVED	0-15 LPM	0	1/2	1	1 1/2	2	3	4	5	6	8	10	15	5/CASE
36353	BARB OUTLET • BRASS-SLEEVED	0-25 LPM	0	1/2	1	2	3	4	5	6	8	10	15	25	5/CASE
36354	MINI • BARB OUTLET • BRASS-SLEEVED	0-4 LPM	0	1/32	1/16	1/8	1/4	1/2	3/4	1	1 1/2	2	3	4	5/CASE
36355	MINI • BARB OUTLET • BRASS-SLEEVED	0-8 LPM	0	1/8	1/2	1	1 1/2	2	2 1/2	3	4	5	6	8	5/CASE
36356	MINI • BARB OUTLET • BRASS-SLEEVED	0-15 LPM	0	1/2	1	1 1/2	2	3	4	5	6	8	10	15	5/CASE
36357	MINI • BARB OUTLET • BRASS-SLEEVED	0-25 LPM	0	1/2	1	2	3	4	5	6	8	10	15	25	5/CASE
36358	EXTENDED • DISS OUTLET • BRASS-SLEEVED	0-8 LPM	0	1/8	1/2	1	1 1/2	2	2 1/2	3	4	5	6	8	5/CASE
36359	EXTENDED • DISS OUTLET • BRASS-SLEEVED	0-15 LPM	0	1/2	1	1 1/2	2	3	4	5	6	8	10	15	5/CASE
36360	BARB OUTLET • ALL BRASS	0-25 LPM	0	1/2	1	2	4	6	8	10	12	15	20	25	5/CASE
36361	BARB OUTLET • 2 DISS OUTLETS • ALL BRASS	0-25 LPM	0	1/2	1	2	4	6	8	10	12	15	20	25	5/CASE



WE CARE LIKE FAMILY™

1.888 DYNAREX

DYNAREX.COM

10 GLENSHAW STREET | ORANGEBURG, NY 10962

Contact your local sales representative or call to find your nearest distributor.



R210823

©2021 Dynarex Corporation