

WE CARE LIKE FAMILY

Intubating Stylettes

Facilitates quick intubation of endotracheal tubes


- Rigid aluminum wire within a flexible plastic sheath helps reduce friction between the stylette and tracheal tube
- Available in four sizes to fit most adult and pediatric tracheal tubes
- For use in hospitals, nursing homes, and EMS
- Not made with natural rubber latex

Compares to Rusch Flexi-Slip™ Stylet*

For our complete line of **EMS Products**, download our EMS Catalog at Dynarex.com



See our website for our Endotracheal Tube Introducers (Bougie) (#4582-4583)



Reorder No.	Description	Quantity
4584	Intubating stylettes, 6FR	20/Cs
4585	Intubating stylettes, 10FR	20/Cs
4586	Intubating stylettes, 12FR	20/Cs
4587	Intubating stylettes, 14FR	20/Cs

*Rusch Flexi-Slip is a trademark of Teleflex Incorporated.



Find us on:



dynarex[®] is your one-stop medical supply company since 1967

EMS PROFESSIONAL PREFERRED PRODUCTS

For more information contact your local Dynarex® distributor sales representative. Or call **888-396-2739** to find your nearest distributor. 10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarex.com