## **ALOESKIN GLOVES**

WE CARE LIKE FAMILY™

## Nitrile Exam Gloves with Aloe Vera

## **FEATURES AND BENEFITS**

- Nitrile exam gloves specially coated with Aloe Vera
- Aloe Vera and Vitamin E nourish and protect skin
- Textured fingertips and beaded cuffs
- Ambidextrous and powder-free
- 2 mil. minimum thickness
- Not made with natural rubber latex
- **HCPCS**: A4927







Reorder No.	Product Name	Quantity
6721	AloeSkin Nitrile Exam Gloves – Small	10/100/Cs
6722	AloeSkin Nitrile Exam Gloves – Medium	10/100/Cs
6723	AloeSkin Nitrile Exam Gloves – Large	10/100/Cs
6724	AloeSkin Nitrile Exam Gloves – X-Large	10/100/Cs















**dynarex**<sup>®</sup> is your one-stop medical products brand since 1967

For more information contact your local Dynarex sales representative OR Call us at: (888) 396-2739 • www.dynarex.com