



WE CARE LIKE FAMILY™



## Sphygmomanometers

**Great for medical professionals to monitor blood pressure.**

- Accurate blood pressure readings with state of the art aneroid sphygmomanometer
- Chrome-plated, precision crafted, easy to read 300 mmHg with air release valve
- Ideal for the cost-conscious healthcare provider looking for quality and affordability
- Not made with natural rubber latex
- **Key and instruction book are contained within a small plastic bag inside the canvas bag**

Reorder No.	Sphygmomanometer	Quantity
7104	Infant (Under 1 year old)	10/Cs
7105	Toddler (1-4 years old)	10/Cs
7106	Child (5-9 years old)	10/Cs
7107	Adult /Medium (Arm)	10/Cs
7108	Adult /Large (Arm)	10/Cs
7109	Adult /XL (Thigh)	10/Cs

Reorder No.	Description	Quantity	Cuff Size	Min. Arm Size	Max. Arm Size
7104	Infant (Under 1 year old)	10/Cs	26x4.75 cm (10x 2")	12 cm (4.7")	19 cm (7.5")
7105	Toddler (1-4 years old)	10/Cs	33x5 cm (13x2")	17 cm (4.7")	26 cm (10.25")
7106	Child (5-9 years old)	10/Cs	38x7 cm (15x2.75")	17 cm (4.7")	27 cm (10.5")
7107	Adult /Medium (Arm)	10/Cs	51x14 cm (20x5.5")	27 cm (4.7")	36 cm (14.25")
7108	Adult /Large (Arm)	10/Cs	62x17.5 cm (25x7")	32 cm (4.7")	46 cm (18")
7109	Adult /XL (Thigh)	10/Cs	72x21 cm (28x8.25")	38 cm (4.7")	53 cm (21")

See our website for our Single and Dual Head Stethoscopes (7110 & 7111)



\*Depending on market demand and Dynarex stock, we will be happy to address any orders for items currently not shown.



Find us on:



**dynarex**® is your one-stop medical supply company since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • [www.dynarex.com](http://www.dynarex.com)