



WE CARE LIKE FAMILY™

## SensiLance™ Pressure & Button Activated Safety Lancets

### For Painless Sampling

#### SensiLance™ Pressure Activated Lancets

- For capillary blood samples
- Ready to use - no loading needed
- Needle automatically retracts
- Not made with natural rubber latex

#### SensiLance™ Button Activated Lancets

- Tri-bevel tip for virtually painless sampling
- Needle automatically retracts
- For capillary blood samples
- No lancing device required
- Not made with natural rubber latex



Reorder No.	Description	Quantity
7112	21 gauge - 2.2 mm Depth, Pressure Activated	10/100/Cs
7114	26 gauge - 1.8 mm Depth, Pressure Activated	10/100/Cs
7115	28 gauge - 1.8 mm Depth, Pressure Activated	10/100/Cs
7122	21 gauge - 2.4 mm Depth, Button Activated	10/100/Cs
7124	26 gauge - 1.8 mm Depth, Button Activated	10/100/Cs
7125	28 gauge - 1.8 mm Depth, Button Activated	10/100/Cs

7112-7125 HCPC: A4259

**Other Products Available:**  
See our website for our complete line of exam gloves



Find us on:



**dynarex®** is your one-stop medical products brand since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • [www.dynarex.com](http://www.dynarex.com)