



WE CARE LIKE FAMILY

Pen Needles

Used in conjunction with injection pens to deliver injectable medications into the body

- Come in a variety of needle lengths and diameters and are used by both health professionals and patients
- Non Toxic, Pyrogen free
- Not made with natural rubber latex
- Sterile
- Available in 30G-5mm, 30G-8mm, 31G-5mm, 31G-8mm, and 32G-4mm

Compares to:
Ultimed,
MHC Medical Products,
BD, Owen Mumford



Reorder No.	Description	Quantity
NEW! 7024	Pen Needles, 30G x 5mm	12/100/Cs
NEW! 7025	Pen Needles, 31G x 8mm	12/100/Cs
7026	Pen Needles, 30G x 8mm	12/100/Cs
7027	Pen Needles, 31G x 5mm	12/100/Cs
7028	Pen Needles, 32G x 4mm	12/100/Cs

Compatible with these insulin pens:

- Autopen®
- BD pen
- FlexPen®
- NovoPen® / NovoLet® / InnoLet®
- HumaPen Ergo®
- HumaPen® Luxura™ /
- HumaPen® Memoir™
- KwikPen™
- Solostar®
- SymlinPen®
- BerliPen®
- Omnican Pen®
- OptiClik™
- OptiPen® / OptiSet®

See our website for our full line of Insulin Syringes (6901-6905, 6940-6959)



Find us on:



dynarex® is your one-stop medical supply company since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarex.com