

Single Post *Stainless Steel Mayo Stand



MODEL: MS-23

- Stainless steel frame and base
- Tip style position base
- Removable stainless steel tray
- 2 ball bearing rubber casters
- 1" pole and frame
- Knob lock height adjustment
- Removable 19" x 12-3/4" x 3/4" stainless steel tray
- 15 lbs. load capacity under normal use

*Clinton stainless steel items are not M.R.I. compatible.

Shipping Weight: 17 lbs.

SPECIFICATIONS

Base: 18" W x 22" D

Height Range: 31"—50"

California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1