



Stainless Steel Medicine Cup 2 oz 4241 12/BX

Product Features	Made of the high grade stainless steel (18/8,304)	
	Crafted without seams which prevents bacteria growth and allows for easy cleaning	
	Autoclavable	

Attributes	Brand Name	Tech-Med
	Size	2 1/4" x 2 1/16"
	Material	18/8 Stainless Steel
	Sterile	No
	Weight	1.30073 lbs
Packaging	Packaging	12/BX
	Case Packaging	12/BX
	Total Case Quantity	12
	Item Packaging Barcode	665973021843
	Inner Packaging Barcode	(01)10665973021840
	Case Packaging	(01)20665973021847

Barcode

Phone: 1 (631) 656-3800 Toll Free: 1 (800) 243-0741 Fax: 1 (631) 656-3810

