



Steel Drop-Arm Commode with Wheels and Padded Armrests

HCPCS: E0165 OR E0166









FIGURE A

FIGURE B

B FIGURE C

URE C FIGURE D









FIGURE E

FIGURE F

FIGURE G

FIGURE H

- Easy-to-release arm mechanism allows for safe lateral patient transfers to and from commode
- Arm provides additional support while transferring
- Manufactured with 1" steel
- White powder coated frame is attractive and easy to maintain
- Padded arm rests provide extra comfort and support
- Comes with 7.5 qt. commode bucket with carry handle and splash shield
- Rust resistant, swivel casters (2 with locks)
- Limited Lifetime Warranty
- Made using recycled materials

| PRODUCTS | Description |
|------------------------|---------------------------------|
| 11101W-2 | Commode,Drop Arm w/ Wheels,2/cs |
| SPECIFICATIONS | Description |
| Carton Shipping Weight | 47 lbs. |
| Height | 17" - 21" |
| Seat Dimensions | 14.25" (W) x 16.75" (D) |
| Warranty | Limited Lifetime |
| Weight | 22 lbs. |
| Weight Capacity | 300 lbs |
| Width | 18" between arms 22" Outside |