



Products / Traditional Wound Care / Gauze Dressings / Super Sponges / Super Sponges Medium 32-Ply



## Super Sponges Medium 32-Ply

2334

100/BG, 6 BG/CS

### Product Features

Multi-layer washed, fluff-dried, 100% cotton

Highly absorbent and durable

Thick, diagonally folded construction adds loft to enhance its overall performance

### Attributes

<b>Brand Name</b>	Dukal
<b>Size</b>	Medium
<b>Ply</b>	32-Ply
<b>HCPC Code</b>	A6216
<b>Sterile</b>	No
<b>Weight</b>	6.1068 lbs

### Packaging

<b>Packaging</b>	100/BG, 6 BG/CS
<b>Inner Packaging</b>	100/BG
<b>Case Packaging</b>	6 BG/CS
<b>Total Case Quantity</b>	600
<b>Inner Packaging Barcode</b>	(01)10665973002788
<b>Case Packaging Barcode</b>	(01)20665973002785

2 Fleetwood Ct.  
Ronkonkoma, NY 11779

**Phone:** 1 (631) 656-3800  
**Toll Free:** 1 (800) 243-0741  
**Fax:** 1 (631) 656-3810

