



## Surgical Mask with Ear Loop 3-Ply, Blue

1541

Ear Loop, 50/BX, 6 BX/CS

### Product Features

3-Ply pleated lightweight, breathable material

BFE of  $\geq 99\%$

Meets ASTM Level 1 standard

Available in earloop or tie on with flexible nosepiece

### Attributes

<b>Brand Name</b>	Dukal
<b>Color</b>	Blue
<b>HCPC Code</b>	A4928
<b>Sterile</b>	No
<b>Weight</b>	3.96832 lbs

### Packaging

<b>Packaging</b>	50/BX, 6 BX/CS
<b>Inner Packaging</b>	50/BX
<b>Case Packaging</b>	6 BX/CS
<b>Total Case Quantity</b>	300
<b>Inner Packaging Barcode</b>	(01)10665973005550
<b>Case Packaging Barcode</b>	(01)20665973005557