

## **Economy Step Stool with Hand Rail**



## **MODEL: T-50**

- All welded steel base and top (no screws or rivets)
- Rubber tread top for safety
- Reinforced rubber feet
- Durable chrome finish
- 7/8" diameter handrail 36" high
- 350 lbs. load capacity under normal use

Shipping weight: 13 lbs.

**SPECIFICATIONS** 

PLATFORM: 14.25" L X 11.25" W X 9"H



California residents please be advised of the following, pursuant to Proposition 65:

WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.