SafeBasics[®] Dry-Back Bibs



CERTIFICATION & STANDARDS

ISO 13485 certified manufacturer. Health Canada registered medical device manufacturer.

STORAGE CONDITIONS

Bibs should be disposed of after each use. Store in a cool, dry place.



S Medicom[®] TECHNICAL DATA SHEET

FEATURES & BENEFITS

2-Ply Bibs

- o Durable, highly absorbent patient protection.
- Horizontal embossed pattern and unique water-repellent edge for enhanced protection from fluids.
- Strong, tear-resistant plastic backing increases clothing protection.
- Suitable for cleanings and minor procedures.
- Available in a variety of colours to coordinate with Medicom SafeBasics® plastic cups.

RECOMMENDED USE

• Bibs are made of paper intended to absorb saliva or moisture from the oral cavity during procedures. Latex-free.

PRODUCT DETAILS

Product Name: Product Type:

Size:

Quantity per Case: Color: Country of Origin: Medicom SafeBasics® Dry-Back Bibs Single-use, non-sterile, 2 tissue ply + 1 PE poly 13" x 18" (33 x 45,7 cm) 500 Bibs See Appendix 1 USA

BIB SPECIFICATIONS						
CHARACTERISTIC	SPECIFICATION					
Material Composition	Tissue Paper & Polyethylene Film					
Description	Folded piece of paper with a plastic backing. Embossed horizontal pattern on the body of the bib and vertical embossed patter on the vertical borders.					
	Unfolded Size	Unfolded Width	Unfolded Length			
Dimension	13" x 18" (33 x 45, 7 cm)	13″	18"			
	Folded Size	Folded Width	Folded Length			
	4 ⁵ / ₈ " x 6 ⁵ / ₈ " (11,8 x 16,8 cm)	4 ⁵ / ₈ "	6 ⁵ / ₈ "			



Appendix 1 (List of items)								
Item Code	Description	Packaging Information	Size	Color	Specific Attribute	Treatment	Marketed in USA as	Marketed in Canada as
8280	Medicom SafeBasics Dry-Back Bibs	500 Bibs/Case	13" x 18" (33 x 45,7 cm)	Lavender	2 Ply	Non-Sterile	General Purpose	General Purpose
8281	Medicom SafeBasics Dry-Back Bibs	500 Bibs/Case	13" x 18" (33 x 45,7 cm)	Silver Grey	2 Ply	Non-Sterile	General Purpose	General Purpose
8282	Medicom SafeBasics Dry-Back Bibs	500 Bibs/Case	13" x 18" (33 x 45,7 cm)	Blue	2 Ply	Non-Sterile	General Purpose	General Purpose
8283	Medicom SafeBasics Dry-Back Bibs	500 Bibs/Case	13" x 18" (33 x 45,7 cm)	White	2 Ply	Non-Sterile	General Purpose	General Purpose
8284	Medicom SafeBasics Dry-Back Bibs	500 Bibs/Case	13" x 18" (33 x 45,7 cm)	Yellow	2 Ply	Non-Sterile	General Purpose	General Purpose
8285	Medicom SafeBasics Dry-Back Bibs	500 Bibs/Case	13" x 18" (33 x 45,7 cm)	Green	2 Ply	Non-Sterile	General Purpose	General Purpose

Appendix 1 (List of items)								
Item Code	Description	Packaging Information	Size	Color	Specific Attribute	Treatment	Marketed in USA as	Marketed in Canada as
8286	Medicom SafeBasics Dry-Back Bibs	500 Bibs/Case	13" x 18" (33 x 45,7 cm)	Dusty Rose	2 Ply	Non-Sterile	General Purpose	General Purpose
8287	Medicom SafeBasics Dry-Back Bibs	500 Bibs/Case	13" x 18" (33 x 45,7 cm)	Beige	2 Ply	Non-Sterile	General Purpose	General Purpose
8288	Medicom SafeBasics Dry-Back Bibs	500 Bibs/Case	13" x 18" (33 x 45,7 cm)	Peach	2 Ply	Non-Sterile	General Purpose	General Purpose
8289	Medicom SafeBasics Dry-Back Bibs	500 Bibs/Case	13" x 18" (33 x 45,7 cm)	Aqua	2 Ply	Non-Sterile	General Purpose	General Purpose

AMD Medicom Inc. | 2555 Chemin de l'Aviation, Montreal, Quebec, Canada H9P 2Z2