



## TECHNICAL DATA SHEET

### AMD-Ritmed Laboratory Coat Items # A8030 to A8051

#### Intended Use

Disposable SMS Laboratory Coat designed to protect health professionals, lab technicians and other users from fluids and cross contamination.

Characteristic	Specification		
Description	Coat made of non-woven material with snap button front closure, knitted collar and cuffs. No pockets.		
Material	Spunbond Meltblown Spunbond (SMS) Polypropylene. Latex-Free		
Packaging	Gowns / Bag	Bags / Case	Gowns / Case
	10	5	50

#### Item Code Table

Color	Size									
	Small	Medium	Large	XL	XXL	3XL	4XL	7XL	9XL	12XL
Bleu	A8031	A8032	A8033	A8034	A8035	A8036	A8038	A8039	A8040	A8030
White	A8041	A8042	A8043	A8044	A8045	A8046	A8048	A8049	A8050	A8051

#### Dimension Table

Dimension											
Measurements	Raglan Sleeve Style						Straight Sleeve Style				
	Small Size (cm)	Medium Size (cm)	Large Size (cm)	XL Size (cm)	XXL Size (cm)	3XL Size (cm)	4XL Size (cm)	7XL Size (cm)	9XL Size (cm)	12XL Size (cm)	
(A) Underarm Width	62,5	64,5	66,5	69,0	72,0	73,5	77,0	86,0	95,0	101,0	
(B) Lower Sleeve Length	48,0	49,0	51,0	52,0	53,0	55,0	57,5	65,0	70,0	68,0	
(C) Cuff Length	5,0	5,0	5,0	5,0	5,0	5,0	5,0	5,0	5,0	5,0	
(D) Shoulder Width	7,0	7,0	7,0	7,0	7,0	7,0	7,0	7,0	7,0	7,0	
(E) Neck Width (stretched)	19,0	19,0	21,0	21,0	23,0	23,0	24,0	25,0	26,0	27,0	
(F) Open Gown Width	134	139	144	149	154	158	164	182	198	212	
(G) Gown Length	93,5	95,0	97,5	103,0	107,5	112,0	114,5	121,5	146,0	146,0	
(H) Shoulder Height*	31,0	32,0	32,0	33,0	35,0	35,0	31,0	34,0	34,0	34,0	
(I) Back Slit	30,5	30,5	30,5	30,5	30,5	30,5	30,5	30,5	30,5	30,5	
(J) Sleeve Width at Wrist (stretched)	15,0	15,0	15,0	15,0	15,0	15,0	15,0	15,0	15,0	15,0	
(K) Cuff Width	5,0	5,0	5,0	5,0	5,0	5,0	5,0	5,0	5,0	5,0	
- Number of Buttons	5	5	5	5	5	5	6	6	6	6	

TECHNICAL DATA SHEET

AMD-Ritmed Laboratory Coat  
Items # A8030 to A8051

Schematic representation of “raglan sleeve” gown (back view)





Figure 2 – Schematic representation of “straight sleeve” gown (back view)





## TECHNICAL DATA SHEET

### AMD-Ritmed Laboratory Coat Items # A8030 to A8051

Item Code	Description	Color	Specific Attributes
A8030	AMD-Ritmed Laboratory Coat 12XL Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	Blue	Non-Sterile, Back Slit, Straight Sleeve
A8031	AMD-Ritmed Laboratory Coat Small Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	Blue	Non-Sterile, Back Slit, Raglan Sleeve
A8032	AMD-Ritmed Laboratory Coat Medium Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	Blue	Non-Sterile, Back Slit, Raglan Sleeve
A8033	AMD-Ritmed Laboratory Coat Large Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	Blue	Non-Sterile, Back Slit, Raglan Sleeve
A8034	AMD-Ritmed Laboratory Coat XL Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	Blue	Non-Sterile, Back Slit, Raglan Sleeve
A8035	AMD-Ritmed Laboratory Coat XXL Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	Bleu	Non-Sterile, Back Slit, Raglan Sleeve
A8036	AMD-Ritmed Laboratory Coat 3XL Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	Blue	Non-Sterile, Back Slit, Raglan Sleeve
A8038	AMD-Ritmed Laboratory Coat 4XL Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	Blue	Non-Sterile, Back Slit, Straight Sleeve
A8039	AMD-Ritmed Laboratory Coat 7XL Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	Blue	Non-Sterile, Back Slit, Straight Sleeve
A8040	AMD-Ritmed Laboratory Coat 9XL Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	Blue	Non-Sterile, Back Slit, Straight Sleeve
A8041	AMD-Ritmed Laboratory Coat Small Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	White	Non-Sterile, Back Slit, Raglan Sleeve
A8042	AMD-Ritmed Laboratory Coat Medium Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	White	Non-Sterile, Back Slit, Raglan Sleeve
A8043	AMD-Ritmed Laboratory Coat Large Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	White	Non-Sterile, Back Slit, Raglan Sleeve
A8044	AMD-Ritmed Laboratory Coat XL Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	White	Non-Sterile, Back Slit, Raglan Sleeve
A8045	AMD-Ritmed Laboratory Coat XXL Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	White	Non-Sterile, Back Slit, Raglan Sleeve
A8046	AMD-Ritmed Laboratory Coat 3XL Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	White	Non-Sterile, Back Slit, Raglan Sleeve
A8048	AMD-Ritmed Laboratory Coat 4XL Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	White	Non-Sterile, Back Slit, Straight Sleeve



## TECHNICAL DATA SHEET

### AMD-Ritmed Laboratory Coat Items # A8030 to A8051

Item Code	Description	Color	Specific Attributes
A8049	AMD-Ritmed Laboratory Coat 7XL Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	White	Non-Sterile, Back Slit, Straight Sleeve
A8050	AMD-Ritmed Laboratory Coat 9XL Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	White	Non-Sterile, Back Slit, Straight Sleeve
A8051	AMD-Ritmed Laboratory Coat 12XL Size, 10 Gowns/Bag, 5 Bags/Case, 50 Gowns/Case	White	Non-Sterile, Back Slit, Straight Sleeve

#### Regulatory

FDA Class 1, FDA Product Code OEA  
Health Canada Class 1

#### Storage Conditions

Store under ambient temperature and dry conditions

#### Country of origin

China