


# TECHNICAL DATA SHEET Plastic Tipped Applicator

#### Intended Use

Medical - General Hospital

Characteristic	Specification						
Material Composition	Shaft	Plastic Products – Hollow Polystyrene					
	Bud	Polyester - Medical Grade, First Quality Polyester Dull Finish					
	Adhesive	Non-Toxic Sterilizable Hot Adhesive					
	Binder	Bud is formed and held, using a non-toxic, water soluble adhesive applied evenly over the bud in a mist application. (1% Single Methocel)					
Description	Serrated shaft with soft fibrous material at one end.						
Color	White						
Dimension	Bud	Length (L2)		Diameter (D1)			
		0.625"		0.218"			
	Swab	Overall Length (L1)		Overhang Length (L3)			
		6"		0.1056"			
	Shaft	Length (L4)	Outer Diameter (D2)		Inner Diameter (D2)		
		5.875"	0.100"		0.060"		



#### Additional Information

Storage Condition	Country of Origin		
Add recommended storage condition	Add country of origin		

#### Appendix 1 List of Item Codes and Attributes

TDS 07.01-13EN - A Date: 06-Aug-2020 Page 1 of 1 FRM 5.04B-1 Rev 03

## Appendix 1 (List of items)

### AMD-Medicom, Inc.

Item Code	Description	Packaging Information	Size	Treatment	Marketed in USA as	Marketed in Canada as
4001D	Purfybr Plastic Tipped Applicator	1 Applicator/Pouch, 100 Pouches/Box, 1000 Applicators/Case; Single wrapped in paper	6"	Sterile	Medical Device Class 1	Medical Device Class 1
4002D	Purfybr Plastic Tipped Applicator	2 Applicators/Pouch, 100 Pouches/Box, 2000 Applicators/Case; Double wrapped in paper	6"	Sterile	Medical Device Class 1	Medical Device Class 1

APP1 07.01-13EN\_A Date: 06-Aug-2020 Page 1 / 1