

## TECHNICAL DATA SHEET Cohesive Bandage Roll

### Intended Use

Cohesive bandage rolls used in non-invasive procedures such as the covering of a dressing or to support a body part as a result of an injury. Used in health care facilities.

Characteristic	Specification		
Material Composition	<b>Latex</b>	Flexible non-woven cloth (60- 70%): 84% Polypropylene / 12.5% Spandex / 3.5% Natural Latex	
		Adhesive Latex (30 - 40%): 91% Natural Latex / 5% White Oil / 2% Antioxidants	
	<b>Latex-Free</b>	Flexible non-woven cloth (75-90%): 69-84% Polypropylene / 6% Spandex / 10-25% Synthetic Glue	
Description	Self-adherent bandage roll		
Dimension	<b>Size</b>	<b>Width (cm)</b>	<b>Length (m)</b>
	1" x 5 yds	2,5	≥ 4,5
	2" x 5 yds	5,1	
	3" x 5 yds	7,6	
	4" x 5 yds	10,2	
	6" x 5 yds	15,3	

### Additional Information

Regulatory	Storage Condition	Country of Origin
FDA Class 1 Health Canada Class 1	Store under ambient temperature and dry condition	China

### Appendix 1 List of Item Codes and Attributes

## Appendix 1 (List of items)

AMD Medicom Inc.

Item Code	Description	Packaging Information	Size	Color	Specific Attribute	Treatment
A1061-A	AMD-Ritmed Cohesive Bandages	72 Rolls/Box, 6 Boxes/Case, 432 Rolls/Case Individually Wrapped	1" x 5 yds (2,5 cm x 4,5 m)	Assorted colors (6 colors), Brigh Blue, Pink, Ligth Purple, Hot Orange, Yellow & Lime Green	Latex	Non-Sterile
A1061-B	AMD-Ritmed Cohesive Bandages	72 Rolls/Box, 6 Boxes/Case, 432 Rolls/Case Individually Wrapped	1" x 5 yds (2,5 cm x 4,5 m)	Blue	Latex	Non-Sterile
A1061-T	AMD-Ritmed Cohesive Bandages	72 Rolls/Box, 6 Boxes/Case, 432 Rolls/Case Individually Wrapped	1" x 5 yds (2,5 cm x 4,5 m)	Tan or Skin	Latex	Non-Sterile
A1061-TLF	AMD-Ritmed Cohesive Bandages	72 Rolls/Box, 6 Boxes/Case, 432 Rolls/Case Individually Wrapped	1" x 5 yds (2,5 cm x 4,5 m)	Tan or Skin	Latex-free	Non-Sterile
A2061-A	AMD-Ritmed Cohesive Bandages	36 Rolls/Box, 6 Boxes/Case, 216 Rolls/Case Individually Wrapped	2" x 5 yds (5 cm x 4,5 m)	Assorted colors (6 colors), Brigh Blue, Pink, Ligth Purple, Hot Orange, Yellow & Lime Green	Latex	Non-Sterile
A2061-ALF	AMD-Ritmed Cohesive Bandages	36 Rolls/Box, 6 Boxes/Case, 216 Rolls/Case Individually Wrapped	2" x 5 yds (5 cm x 4,5 m)	Assorted colors (6 colors), Brigh Blue, Red, Green, Pink, Ligth Purple, & White	Latex-free	Non-Sterile

Item Code	Description	Packaging Information	Size	Color	Specific Attribute	Treatment
A2061-B	AMD-Ritmed Cohesive Bandages	36 Rolls/Box, 6 Boxes/Case, 216 Rolls/Case Individually Wrapped	2" x 5 yds (5 cm x 4,5 m)	Blue	Latex	Non-Sterile
A2061-MIX	AMD-Ritmed Cohesive Bandages	36 Rolls/Box, 6 Boxes/Case, 216 Rolls/Case Individually Wrapped	2" x 5 yds (5 cm x 4,5 m)	Assorted designs: 8 Bones, 7 Hearts, 7 Stars, 7 Blue camouflage, 7 Purple camouflage	Latex	Non-Sterile
A2061-PPA	AMD-Ritmed Cohesive Bandages	36 Rolls/Box, 6 Boxes/Case, 216 Rolls/Case Individually Wrapped	2" x 5 yds (5 cm x 4,5 m)	Paw printed assorted colors	Latex	Non-Sterile
A2061-R	AMD-Ritmed Cohesive Bandages	36 Rolls/Box, 6 Boxes/Case, 216 Rolls/Case Individually Wrapped	2" x 5 yds (5 cm x 4,5 m)	Red	Latex	Non-Sterile
A2061-T	AMD-Ritmed Cohesive Bandages	36 Rolls/Box, 6 Boxes/Case, 216 Rolls/Case Individually Wrapped	2" x 5 yds (5 cm x 4,5 m)	Tan or Skin	Latex	Non-Sterile
A2061-TLF	AMD-Ritmed Cohesive Bandages	36 Rolls/Box, 6 Boxes/Case, 216 Rolls/Case Individually Wrapped	2" x 5 yds (5 cm x 4,5 m)	Tan or Skin	Latex-free	Non-Sterile
A3061-A	AMD-Ritmed Cohesive Bandages	24 Rolls/Box, 6 Boxes/Case, 144 Rolls/Case Individually Wrapped	3" x 5 yds (7,5 cm x 4,5 m)	Assorted colors (6 colors), Brigh Blue, Pink, Ligth Purple, Hot Orange, Yellow & Lime Green	Latex	Non-Sterile

Item Code	Description	Packaging Information	Size	Color	Specific Attribute	Treatment
A3061-ALF	AMD-Ritmed Cohesive Bandages	24 Rolls/Box, 6 Boxes/Case, 144 Rolls/Case Individually Wrapped	3" x 5 yds (7,5 cm x 4,5 m)	Assorted colors (6 colors), Brigh Blue, Red, Green, Pink, Ligth Purple, & White	Latex-free	Non-Sterile
A3061-T	AMD-Ritmed Cohesive Bandages	24 Rolls/Box, 6 Boxes/Case, 144 Rolls/Case Individually Wrapped	3" x 5 yds (7,5 cm x 4,5 m)	Tan or Skin	Latex	Non-Sterile
A3061-TLF	AMD-Ritmed Cohesive Bandages	24 Rolls/Box, 6 Boxes/Case, 144 Rolls/Case Individually Wrapped	3" x 5 yds (7,5 cm x 4,5 m)	Tan or Skin	Latex-free	Non-Sterile
A4061-A	AMD-Ritmed Cohesive Bandages	18 Rolls/Box, 6 Boxes/Case, 108 Rolls/Case Individually Wrapped	4" x 5 yds (10 cm x 4,5 m)	Assorted colors (6 colors), Brigh Blue, Pink, Ligth Purple, Hot Orange, Yellow & Lime Green	Latex	Non-Sterile
A4061-ALF	AMD-Ritmed Cohesive Bandages	18 Rolls/Box, 6 Boxes/Case, 108 Rolls/Case Individually Wrapped	4" x 5 yds (10 cm x 4,5 m)	Assorted colors (6 colors), Brigh Blue, Red, Green, Pink, Ligth Purple, & White	Latex-free	Non-Sterile
A4061-B	AMD-Ritmed Cohesive Bandages	18 Rolls/Box, 6 Boxes/Case, 108 Rolls/Case Individually Wrapped	4" x 5 yds (10 cm x 4,5 m)	Blue	Latex	Non-Sterile
A4061-PPA	AMD-Ritmed Cohesive Bandages	18 Rolls/Box, 6 Boxes/Case, 108 Rolls/Case Individually Wrapped	4" x 5 yds (10 cm x 4,5 m)	Paw printed assorted colors	Latex	Non-Sterile

Item Code	Description	Packaging Information	Size	Color	Specific Attribute	Treatment
A4061-R	AMD-Ritmed Cohesive Bandages	18 Rolls/Box, 6 Boxes/Case, 108 Rolls/Case Individually Wrapped	4" x 5 yds (10 cm x 4,5 m)	Red	Latex	Non-Sterile
A4061-T	AMD-Ritmed Cohesive Bandages	18 Rolls/Box, 6 Boxes/Case, 108 Rolls/Case Individually Wrapped	4" x 5 yds (10 cm x 4,5 m)	Tan or Skin	Latex	Non-Sterile
A4061-TLF	AMD-Ritmed Cohesive Bandages	18 Rolls/Box, 6 Boxes/Case, 108 Rolls/Case Individually Wrapped	4" x 5 yds (10 cm x 4,5 m)	Tan or Skin	Latex-free	Non-Sterile
A4061-TLFB	AMD-Ritmed Cohesive Bandages	18 Rolls/Box, 6 Boxes/Case, 100 Rolls/Case Individually Wrapped	4" x 5 yds (10 cm x 4,5 m)	Tan or Skin	Latex-free	Non-Sterile
A6061-A	AMD-Ritmed Cohesive Bandages	12 Rolls/Box, 6 Boxes/Case, 72 Rolls/Case Individually Wrapped	6" x 5 yds (15 cm x 4,5 m),	Assorted colors (6 colors), Brigh Blue, Pink, Ligth Purple, Hot Orange, Yellow & Lime Green	Latex	Non-Sterile
A6061-T	AMD-Ritmed Cohesive Bandages	12 Rolls/Box, 6 Boxes/Case, 72 Rolls/Case Individually Wrapped	6" x 5 yds (15 cm x 4,5 m)	Tan or Skin	Latex	Non-Sterile
A6061-TLF	AMD-Ritmed Cohesive Bandages	72 Rolls/Case Individually Wrapped	6" x 5 yds (15 cm x 4,5 m)	Tan or Skin	Latex-free	Non-Sterile
A6061-TLFB	AMD-Ritmed Cohesive Bandages	12 Rolls/Box, 6 Boxes/Case, 72 Rolls/Case Individually Wrapped	6" x 5 yds (15 cm x 4,5 m)	Tan or Skin	Latex-free	Non-Sterile

Item Code	Description	Packaging Information	Size	Color	Specific Attribute	Treatment
404597	AMD-Ritmed Cohesive Bandages	100 Rolls/Case Individually Wrapped	4" x 5 yds (10 cm x 4,5 m)	Lime Green	Latex-free	Non-Sterile