

TECHNICAL DATA SHEET
Krimped Gauze Bandage Roll

Intended Use

Crimped gauze bandage is used to cover and protect wounds from infection and further injuries. This type of product fits all parts of the body without restrictive movement. To be used in health care facilities.

Characteristic	Specification	
Material Composition	100% Cotton, Latex free	
Description	Crimped gauze roll (Refer to Appendix 1 specific attribute)	
Dimension	Width (cm)	Stretched Length (cm)
	11,4	≥ 373 (3.8 yd)

Additional Information

Storage Condition	Country of Origin
Store at ambient temperature and dry conditions	China

Appendix 1 List of Item Codes and Attributes

Appendix 1 (List of items)

AMD Medicom Inc.

Item Code	Description	Packaging Information	Size	Color	Specific Attribute	Treatment
A454	AMD-Ritmed Krimped Gauze Bandage Roll	100 Rolls/Bag, 1 Bag/Case, 100 Rolls/Case Bulk	4.5" x 4.1 yd (11 cm x 3,7 m)	White	6 Ply	Non-Sterile
A455	AMD-Ritmed Krimped Gauze Bandage Roll	1 Roll/Pouch, 100 Pouches/Case, 100 Rolls/Case	4.5" x 4.1 yd (11 cm x 3,7 m)	White	6 Ply	Sterile