

TECHNICAL DATA SHEET

Transparent Dressing

Intended Use

Sterile transparent adhesive dressings consisting of a plastic strip coated on one side with an adhesive. This device is used in health care facilities to protect wounds from the introduction of bacteria and contaminants, in the treatment of skin tears, cuts and abrasions, and to secure object to the skin such as IV lines. Conforms to areas that are often difficult to dress, allows observation of the wound at all times, water and fluid resistant. Latex free.

Characteristic	Specification				
	Backing Plastic Film	Polyethylene			
Material Composition	Substrate	Transparent Polyurethane			
	Adhesive	Acrylic glue			
Description	Die-cut adhesive plastic film covered by two paper strips on the adhesive side.				
Dimension	Size	Width (mm)	Length (mm)		
	4" x 4 3/4" (10 X 12 cm)	100	120		

Additional Information

Storage Condition	Country of Origin		
Avoid extreme temperatures (maintain between – 5°C and +40°C). Keep dry and away from sunlight.	China		

Date: 25-Feb-2021

Appendix 1 List of Item Codes and Attributes

Appendix 1 (List of items)

AMD Medicom Inc.

Item Code	Description	Packaging Information	Size	Color			Marketed in Canada as
A60044	AMD-Ritmed	1 Dressing/Pouch, 50 Pouches/Box, 4	4"x 4 ¾" (10 x 12 cm)	Transparent	Sterile	Medical	Medical
	Transparent	Boxes/Case, 200 Dressings/Case				Device	Device
	Dressing					Class 1	Class 1

APP1 02.10-09EN_1A Date: 25-Feb-2021 Page 1 / 1