



TECHNICAL DATA SHEET Dental Chair Light Shield Barrier

Intended Use

Individual disposable light shield barrier designed exclusively for A-dec LED dental chair light. Medical device designed to reduce cross contamination, speed treatment room and increase the patient's perception of cleanliness.

Characteristic	Specification		
Material Composition	Film	Polyester	
	Adhesive Liner	Polyester	
	Adhesive	UV cured acrylic adhesive	
Description	Transparent barrier with adhesive protected with backing on the four sides locations (Refer to Appendix 1 for specific attribute)		
Color	Clear		
Dimension	Part	Width	Length
	Overall Shield	6.637" (168,6 mm)	7.875" (200 mm)
	Top of Vertical Adhesive	0.562" (14,3 mm)	3.5" (88,9 mm)
	Bottom of Vertical Adhesive		3.0" (76,2 mm)
	Top of Side Adhesive	0.4375" (11,1 mm)	2.75" (69,9 mm)
	Bottom of Side Adhesive		2.25" (57,2 mm)
Position of adhesive	Side	Distance from Center of Barrier Film	
	Top and bottom	3.25" (82,6 mm)	
	Side	2.881" (73,2 mm)	

Additional Information

Storage Condition	Country of Origin
Store in a cool dry place	USA

Appendix 1 List of Item Codes and Attributes

Appendix 1 (List of items)

AMD Medicom Inc.

Item Code	Description	Packaging Information	Specific Attribute	Treatment
9565	Medicom SafeShield Disposable Protective Barriers Designed Exclusively for the A-dec LED Dental Light	10 Shields/Sleeve, 10 Sleeves/Case, 100 Shields/Case	Medicom logo in each corner	Non-Sterile