

SPECEFICATIONS

TILT TABLE

// Heavy Duty Motorized Tilt Table




MODEL	WIDTH	LENGTH	HEIGHT	TILT RANGE
TG2726	31" (78.7 cm)	81.5" (207 cm)	34.5" (87.6 cm)	0 to 90 degrees

FEATURES
500 lb (226.80 Kg) Weight Capacity
UL Recognized Components
Casters – Heavy Duty 5" (12.7 cm) Swivel with Locking Brakes
Wide Footboard for Patient Comfort
Safety Belts - Set of 2, Velcro Closure
Angle Indicator
Accommodates Patient Lifts
Foot Control
2" (5.1 cm) Upholstered Top
Steel Frame

SPECIFICATIONS	
Base:	3" (7.6 cm) Tubular Steel
Finish:	Pearl Grey Standard; Porcelain Optional.
Upholstery:	Imperial Blue, Regimental Blue, Wedgewood, Blue Ridge, Dove, Graphite, Black, Burgundy, Forest, Grotto, Sand, Camel.
Electrical:	120V, 1.6A, 192W, hospital grade plug.
Warranty:	Ten (10) Year Limited Warranty.
*Nominal Dimensions Stated	

POWER SOURCE	NUMBER
VCU	1 Each
Linear Actuator	1 Each

STANDARDS
EN 60601-1/UL 60601-1; On Components
 UL Recognized Components
