

Prevent Falls & Wandering

Early-Warning Chair Belts Know before resident gets up!



For Use With Wheelchair Exit Alarms

These early warning seat belts are designed for use with wheelchairs and all of our Fall Monitors.

Early Warning

The Fall Prevention Wheelchair Seat Belts reduce falls by triggering an alert when the Easy Release Buckle is unbuckled or the Quick Release Hook & Loop strap is opened. This notifies a caregiver before a resident attempts to leave their chair.

Non-Restraining

Both the Easy Release and Quick Release seat belts preserve resident's freedom of movement by not acting as a restraint. The Quick Release Seat Belt is ideal for residents with low dexterity as the Hook & Loop strap can be easily unfastened with one hand.

IMPORTANT: Seat Belts are not designed to restrain or hold individuals in position in their chairs.



Easy Release TL-2109
Adjustable 20" x 54"



Quick Release TL-2109V
Adjustable 20" x 44"

**Custom lengths available*

Key Features:

- Easy Release Seat Belt features a large, red button that is easy to locate and release with one finger
- Quick Release Seat Belt features an easily grasped Hook & Loop strap which unfastens quickly even for those with limited dexterity
- Easy Installation



Call Us Today! (800) 650-3637 | www.smartcaregiver.com