

SPECEIFICATIONS

PARALLEL BARS

// Standard Parallel Bars
Motorized Height & Width Adj.




MODEL	LENGTH
TWG1000.SD	6 ft (1.83 m)
TWG1001.SD	8 ft (2.44 m)
TWG1002.SD	10 ft (3.05 m)
TWG1004.SD	12 ft (3.66 m)
TWG1006.SD	16 ft (4.88 m)
TWG1008.SD	18 ft (5.49 m)

FEATURES
600 lb (272.73 Kg) Weight Capacity
Height & Width Adjustable
Fingertip Controls, Individual or Simultaneous Operation
Steel Frame

SPECIFICATIONS	
Bar Range:	Vertical: 25.5" – 39.75" (64.1 cm - 101 cm) to the top of the bar Horizontal: 15.5" – 24.5" (39.4 cm- 62.2 cm) inside dimension
Width:	Platform: 47 3/16" (119.8 cm) Between Standards: 31.5" (80 cm) Overall Dimensions: 53.5" (135.9 cm)
Base:	10 Gauge Steel; 1 3/4" (4.4 cm) <i>Optional Accessories:</i> 12" (30.5 cm) Tapered End and/or 24" (61 cm) Platform Extension
Handrails:	Anodized Aluminum; 1 5/8" (4.4 cm) Standard O.D.
Finish:	Pearl Grey Standard; Porcelain Optional. Chrome plated Supports.
Electrical:	115V, 6.8A, 782W (230 V/AC, 50 Hz, 3.4A) hospital grade plug
Matting:	Non-skid black mat: width & length to fit.
Warranty:	Ten (10) Year Limited Warranty.
*Nominal Dimensions Stated	

POWER SOURCE	NUMBER
Motors	2 Each

STANDARDS
UL Classified IEC60601-1(ed.3)
 CLASSIFIED c UL US 60601-1
ISO 14971, RM
