



Blood Draw Chair Straight Arm, Blue

4382

1/EA

Product Features



WARNING: This product can expose you to chemicals including chromium (hexavalent compounds), which are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

This quality chair features tool free assembly and padded seat for ultimate patient comfort

One control knob allows you to raise and lower the arms effortlessly (28"-35")

7" wide arm allows for patient comfort and greater flexibility for the phlebotomist

Arms can be raised and lowered in infinite positions and rotate left or right

Welded solid steel frame with a powder coated finish

Vinyl is mildew resistant, sulfide stain resistant and static dissipative

Adjustable levelers on all four legs allow the chair to sit solidly on an uneven floor

500 lb Weight Capacity

Dimensions: Seat Height to Floor 19", Inside Arm Width 21", Seat Depth 17"

Overall dimensions: 38" height, 34" width, 28" depth

Attributes

| | |
|-------------------|--------------|
| Brand Name | Tech-Med |
| Color | Blue |
| Sterile | No |
| Weight | 65.69775 lbs |

Packaging

| | |
|----------------------------|--------------------|
| Packaging | 1/EA |
| Case Packaging | 1/EA |
| Total Case Quantity | 1 |
| Case Packaging | (01)10665973022328 |