

Value That Sticks

# SwiftSet™

## Topical Skin Adhesive

From the Covidien brand you trust with the clinical support you demand.

### Fast

**30-second dry time, 1 layer application**

(strength equivalent to a 4.0 suture<sup>1</sup>)

### Precise

**Two clog-free tip options for pinpoint and wider application**

### Easy

**One-handed activation, flow control, no refrigeration necessary**

*One product for all your clinical needs.*

	SwiftSet™	DermaBond™* Advanced	DermaBond™* Mini
Set Time	30 seconds. Full apposition strength in 1 minute. <sup>2</sup>	Maintain manual approximation for 60 seconds. Full apposition strength in minutes. <sup>3</sup>	1 minute for 2 required applications, with full apposition strength in 2.5 minutes. <sup>4</sup>
Applications Required	1 <sup>2</sup>	1 <sup>3</sup>	2 <sup>4</sup>
Adhesive per Vial	0.8 ml	0.7 ml	0.36 ml
Precision Tip	Yes	No	No
Clog-Free Tip	Yes <sup>2</sup>	No <sup>3</sup>	No <sup>4</sup>
Exothermic Reaction	Low <sup>5</sup>	High <sup>5</sup>	High <sup>5</sup>
Anti-Microbial Barrier	Yes <sup>2</sup>	Yes <sup>3</sup>	Yes <sup>4</sup>
Patient Can Shower Immediately	Yes <sup>2</sup>	Yes <sup>3</sup>	Yes <sup>4</sup>
Violet Tint	Yes	Yes	Yes
Ambient Storage	Yes	Yes	Yes

# SwiftSet™

## Topical Skin Adhesive

Value That Sticks

From the Covidien brand you trust with the clinical support you demand.

### Ordering Information

Reorder Code	Description	Units Per Box
SS6	SwiftSet™ Topical Skin Adhesive	6



### REFERENCES

1. Report on Comparing n-butyl cyanoacrylate skin adhesive (US) Wound Closure Strength Against Competitors. M059 VR 180
2. SwiftSet™ Topical Skin Adhesive Product IFU
3. Dermabond Advanced™ Product IFU
4. Dermabond Product IFU
5. A report on measuring the heat of polymerization of n-butyl cyanoacrylate skin adhesive monomer. Doc No. MLGL-08019-R1

COVIDIEN, COVIDIEN with logo, Covidien logo and *positive results for life* are U.S. and internationally registered trademarks of Covidien AG. All other brands are trademarks of a Covidien company.

™ Trademark of its respective owner.

© 2013 Covidien. 12.13 P130062



555 LONG WHARF DRIVE  
NEW HAVEN, CT  
06511

800-722-8772

WWW.COVIDIEN.COM