

Style Line, Treatment Table with Stool



MODEL 9133

- Features 3 drawers, storage shelf, adjustable backrest and matching stool that stores neatly under the shelf
- Stool is 17.75" wide x 17.75" deep (sits on floor) (45 cm x 45 cm)
- Easy-clean laminate base
- Self-closing drawers with metal sides and 75 lbs. capacity (34 kg)
- Seamless rounded corner top
- 2" firm foam padding
- Heavyweight, knit-backed upholstery
- Also available with non-adjustable backrest (Option 088)
- 400 lbs. load capacity under normal use (181 kg)

SPECIFICATIONS

9133-27	72" L x 27" W x 31" H
9133-30	72" L x 30" W x 31" H

California residents please be advised of the following, pursuant to Proposition 65:

! **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

Options



Laminate Colors



Dark Cherry



Maple



Gray

Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black