



WE CARE LIKE FAMILY™

UNDERCAST PADDING

COTTON

Item # 3880: 2" x 4 yd.

Item # 3881: 3" x 4 yd.

Item # 3882: 4" x 4 yd.

Item # 3883: 6" x 4 yd.

POLYESTER

Item # 3885: 2" x 4 yd.

Item # 3886: 3" x 4 yd.

Item # 3887: 4" x 4 yd.

Item # 3888: 6" x 4 yd.

Dynarex Undercast Padding is designed to provide a cushioning layer under casts for a more comfortable, consistent fit. The Undercast Padding material has a light stretch, making it easy to conform around body contours.

- Designed to provide consistent cushioning between skin and cast
- Crafted from soft, breathable materials for enhanced patient comfort
- Material has a light stretch for easy contouring
- Available in choice of 100% Cotton or Polyester
- Ideal for use with Dynarex Casting Tape
- 20 rolls per bag, 4 bags per case

