


# Venti TENS™ Digital Pain Relief System

The Venti TENS Digital Pain Relief System is a dual-channel prescription digital TENS device equipped with three waveforms and twenty modes of operation to choose from – twelve presets and eight manual. This incredibly advanced, pocket-size portable device is known for its body-part specific settings, which provide ease of application and set-up by patients. Both, its pulse width and pulse frequency are also adjustable.

Used for symptomatic relief and management of chronic pain and acute post-surgical or post-trauma pain.



## FEATURES



12 Body PRESET Programs Include:

- P0- Neck Pain
- P1- Torticollis
- P2- Low Back Pain
- P3- Sciatica
- P4- Epicondylitis
- P5- Wrist Pain
- P6- Post Surgery Knee Pain
- P7- Ankle sprain
- P8- Knee Pain (Gonalgia)
- P9- Hip Pain (Coxalgia)
- P10- Menses Pain (Menstrual)
- P11- Neuralgia

- Dual channel digital \*TENS unit
- Large LCD Display with body-part diagram for targeted pain relief
- Three selectable waveforms: symmetrical biphasic, asymmetrical biphasic and monophasic
- 20 treatment modes- 12 preset and 8 manual
- Last setting recall and settings lock
- Treatment timer and compliance meter
- Protective cover over controls
- Low battery indicator
- Electrode placement suggestions
- Includes: Digital TENS device, pair of lead wires, one package of 4 reusable self-adhering electrodes, 2 AA alkaline batteries, zipper carrying case and manual
- Available as each or in case pack of 20 units

## TECHNICAL SPECIFICATIONS

Channels: Dual channels, isolated between channels

Waveform: Symmetrical Biphasic Rectangular  
Asymmetrical Biphasic Rectangular  
Monophasic Rectangular

Pulse Amplitude: 0 - 130 mA (peak value)

Pulse Width: 50 - 250  $\mu$ s, Adjustable

Pulse Frequency: 1 - 150 Hz, Adjustable

Treatment Time: Continuous, 15, 30, 45, 60, or 90 mins

Dimensions: 6 cm x 2.5 cm x 10 cm

Weight: 5 oz.

Warranty: One year warranty

UPC: 092237615959

## SHIPPING SPECIFICATIONS

Individual Ship Dims: 15" x 10.5" mailer

Individual Ship Weight: 2 lbs.

## CASE PACK SHIPPING SPECIFICATIONS

Case Pack Quantity: 20 units

Case Ship Dims: 14" x 13" x 9" (h x w x d)

Case Ship Weight: 19 lbs.

Item #	Description	Qty
DT2020	Venti TENS Digital Pain Relief System	each

\* Prescription Required: This is a Class II Medical Device. Federal law restricts these devices to sale by or on the order of a licensed health care professional.



**To Order:**  
(800) 3-ROSCOE  
fax: (440) 572-4261  
[roscoemedical.com](http://roscoemedical.com)