


## **VENTURI MASKS**

Dynarex Resp-O₂™ Venturi Masks are designed to deliver oxygen at fixed concentrations ranging from 24% to 50%. These masks come with light green oxygen supply tubing, designed to help improve tube visibility.

- Transparent, molded vinyl mask features smooth edges
- Mask features an adjustable nose clip and elastic strap
- Includes 7' of crush-resistant, star lumen oxygen supply tubing
- Not made with natural rubber latex or DEHP







**Includes** 7' (2.1 m) Crush Resistant **Tubing** 

star Lumen





Select A Vent **Elongated Mask** with 6 Dilutors

STANDARD CONNECTOR

REORDER NO. 33070

ADULT

PACKAGE QUANTITY 50/Case

Select A Vent **Elongated Mask** with 6 Dilutors

STANDARD CONNECTOR

**REORDER NO. 33071** 

PEDIATRIC

PACKAGE QUANTITY 50/Case

Adjustable Dial **Elongated Mask** 

STANDARD CONNECTOR

REORDER NO. 33072

**ADULT** 

PACKAGE QUANTITY 50/Case

Adjustable Dial Elongated Mask


STANDARD CONNECTOR

REORDER NO. 33073

PEDIATRIC

PACKAGE QUANTITY 50/Case





Contact your local







WE CARE LIKE FAMILY