



Retail Packaged Aluminum Crutches with Accessories

HCPCS: E0114

- Lighter than wood, the standard aluminum crutch provides stability and durability
- Double extruded center tube provides additional strength to weight-bearing area
- Pushpin adjustment makes it easy-to-adjust crutch height in 1" increments
- Non-skid, jumbo-size vinyl contoured tips provide excellent traction
- Easy wing nut adjustments for hand grips
- Comfortable, durable accessories are latex free
- Weight Capacity 350 lbs on aluminum crutch (excluding pediatric)

PRODUCTS	Description	
RTL10400	Crutch, Adult,1PR,8/CS	
RTL10402	Crutch,Tall,1PR,8/CS	
SPECIFICATIONS	Description	
Warranty	Limited Lifetime	
Weight Capacity	350 lbs	