



EarthLite®
WORLD'S #1 BRAND IN MASSAGE™

YOSEMITE

TREATMENT TABLE COLLECTION



Yosemite 30



Yosemite, 25"

HEALING HUMANITY THROUGH TOUCH™

YOSEMITE

TREATMENT TABLE COLLECTION

The Yosemite™ massage tables are unlike any Shiatsu massage tables on the market today. The 25" Yosemite's narrow design allows for total client access – perfect for deep tissue bodywork and sports massage. The Yosemite 30's generous 30" width and low height range enables a wide variety of wellness and medical treatments from massage and physical therapy to shiatsu - and ships and stores flat for ease of transport. The powder-coated aluminum adjustable height legs come with Earthlite's patent pending Shimmy-Lok™ leg braces which provide unmatched stability. Topped with a firm and responsive multi-layer cushioning system and wrapped in our supple yet durable Nature's Touch upholstery, the Yosemite is perfect for massage and body treatment therapies. Backed by the industry's best warranty from Earthlite, the #1 brand in Massage since 1987.

Features


- Inviting upholstered round corner top – easy on the therapist, the client and the eyes
- Attractive powder-coated metallic silver legs with black knobs and feet
- Solid wood outlets
- Responsive multi-layer cushioning system (3" on Yosemite 30, 2.5" on Yosemite, 25")
- Round facehole (Yosemite, 25")
- Deluxe Adjustable Facecradle and FacePillow (Yosemite, 30")
- Solid hardwood outlets
- Supple, yet durable, Nature's Touch™ 100% PU upholstery in Black, Agate, or Beige (Beige for Yosemite 30 only)
- Patented Shimmy-Lok™ leg braces

Lifetime Warranty

- Lifetime warranty on the frame and structure
- 3 years on cushioning and upholstery
- CE Mark




NATURE'S TOUCH™ UPHOLSTERY



Compact and Stowable Tables

- Ship flat for cost effective transport and easy storage
- Two easy to install bolts per leg to complete installation (wrench included)



Yosemite Specifications

- Width: 25"
- Length: 74"
- Height Range: 19.25" - 27"
- Static Weight: 3,000 lbs.
- Max Working Weight: 600 lbs.



Yosemite 30 Specifications

- Width: 30"
- Length: 73"
- Height Range: 20" - 26.5"
- Static Weight: 3,000 lbs.
- Max Working Weight: 600 lbs.