



THE LARGEST VOLUME EVACUATED DRAINAGE BOTTLE IS NOW AVAILABLE

IT'S A **BIG** DEAL



BRINGING SAFETY AND CONVENIENCE TO FLUID DRAINAGE MANAGEMENT

ACCEL® Evacuated Drainage Bottles

600mL, 1,000mL and 2,000mL Clear, Shatterproof, Plastic

INTENDED FOR USE TO ASPIRATE, REMOVE, OR SAMPLE BODY FLUIDS

The ACCEL® 600mL, 1,000mL and 2,000mL EVACUATED Drainage Bottles

- Intended for Use to Aspirate, Remove, or Sample Body Fluids
- Safe and Secure, Luer Lock, Needle-Free Access
- Not Made with DEHP, PVC, or Natural Rubber Latex
- Shatterproof, Light Weight, Clear Plastic Bottle Provides Containment of Bio-Hazardous Fluids
- Vacuum up to 26inHq (Inches of Mercury) of Negative Pressure

600mL, 1,000mL and 2,000mL Drainage Bottles



ІТЕМ	Reference Number	Case Qty.
600mL Evacuated Drainage Bottle Individually Packaged, Sterile Fluid Pathway Only	622299	25
1,000mL Evacuated Drainage Bottle without Individual Packaging Sterile Fluid Pathway Only	622275	10
2,000mL Evacuated Drainage Bottle without Individual Packaging Sterile Fluid Pathway Only	622298	6

ACCEL® 600mL, 1,000mL and 2,000mL EVACUATED Drainage Bottles are intended for use to aspirate, remove, or sample body fluids.








Refer to Instructions for Use for relevant warnings, precautions, complications, and contraindications.

Distributed by:

B. Braun Interventional Systems Inc. | Part of the B. Braun Group of Companies 824 Twelfth Avenue | Bethlehem, PA 18018 | USA Tel 877-836-2228 | Fax 610-849-1334 | www.bisusa.org



- Split Closed Male Luer Set
- Spin-Lock® Male Luer Set
- Closed Male Luer Set

Kink-Resistant Tubing for the Aspiration, Removal, or Sampling of Body Fluids Using an External Vacuum Source.

Provides a Safe and Secure Luer Lock Connection Containing Biohazardous Fluids During Drainage. Closed Male Luer Sets Provide for a Safe and Secure Closed System During Evacuated Bottle Exchanges.





ITEM	Reference Number	Length	Case Pack
SPIN-LOCK® SET: Two Spin-Lock® Male Luer Connectors with Pinch Clamp and Non-Vented Spike	470190	51 in. (130cm) approx.	50
CLOSED MALE SET: Closed Male Luer Lock Valve to Spin-Lock® Male Luer Connector with Pinch Clamp	470191	51 in. (130cm) approx.	50
SPLIT CLOSED MALE SET: Bifurcated Set with Two Closed Male Luer Lock Valves to Spin-Lock® Male Luer Connector with Pinch Clamp	470192	48 in. (122cm) approx.	50

ACCEL® Connection Sets are intended for use for the aspiration, removal, or sampling of body fluids using an external vacuum source.





