



Durable Medical Equipment

Moving care forward

Engineer tested, clinician approved.

Because you want to trust the products you use to help care for your patients, we built our durable medical equipment to help you offer quality care. Not only is every product 'clinician-friendly,' they have also been rigorously tested against the benchmarks that matter most, including biocompatibility, weight capacity, fully tested per ISO standards and more.*

Our products are created for clinicians who use mobility products to support their patients every day. We know because they helped us design them. We've standardized our entire DME portfolio, relying heavily on clinical input to ensure we have exactly the right product mix. So you'll spend less time searching—and more time taking care of patients.

This is Cardinal Health today.

And together with you, we're moving care forward.



Asking the experts

We began by working with nurses and physical therapists in the acute, ambulatory and home health markets. Our question was simple: What can we do to improve the products you need, particularly for bariatric patients? Their ideas were practical and powerful. You'll see the results in every product we make.

Here's just a sample:

Need bigger and better wheelchairs to accommodate your patient's needs?

Our wheelchairs offer more than 30 configuration options. No one knows your patients better than you. Configure your wheelchair to fit them best.

Need product labels that are easy to locate and read?

We've moved them where you can see them, so they're highly visible. Our wheelchair labels are clearly visible on the back, instead of hidden underneath the seat.

Saving money two ways

Not only do we offer the right products, our prices are right for the budget-conscious. We also help you control costs by minimizing labor by offering a standardized product set that is quick and easy to manage. What will you do with the time you'll save?

Other ways you'll save time

It's the simple things that really add up. Cardinal Health™ DME products have easy-to-read labels and come with complete product instructions that include tips for safety and cleaning.

All you need, all in one place

Whether you need a new wheelchair, or just spare parts, you'll find the same attention to cost savings and clinical performance.

Our products include:

- **Wheelchairs and accessories**
- **Transport chairs**
- **Walkers and rollators**
- **Crutches and replacement parts**
- **Bathroom accessories**
- **Canes, replacement parts and daily living aids**

When it's time to order, count on our distribution experts to deliver what you need, when you need it. We give you essential tools to do more of what matters most: take care of patients.

*Data on file

An assortment of products to meet patient and caregiver needs.

Wheelchairs and accessories

Take us for a test drive, and you'll never look back.

Because together, we're moving care forward.

Now you get the wheelchairs that fit your

patients best. Need a wide, deep chair with a high weight capacity? Ours offer the sizes and features that are essential to most patients' needs. Want to add an anti-tip device? That's easy. Need a swing-away footrest or desk-length armrest? Just ask.

With over 30 combinations available, we have the wheelchair that's right for your patients. And every one comes standard with the features that matter to clinicians like you—because we asked. The wheelchairs have high quality craftsmanship. The products have been tested extensively. Weight capacity and material number labels are clearly visible and easy to read. All wheelchairs include a carry pocket on the back. The black vinyl upholstery is easy to wipe clean.

Wheelchairs

Cat. no.	Weight capacity	Dimensions (width x depth)	Armrest length type (all removable)	Total width (including wheels)	Packaging
Swing away footrests — chrome finish					
CW0070CS	300 lb.	16 x 16 in.	Desk	25 in.	1 ea/cs
CW0071CS	300 lb.	18 x 16 in.	Desk	27 in.	1 ea/cs
CW0071CSF	300 lb.	18 x 16 in.	Full length	27 in.	1 ea/cs
CW0005CS	300 lb.	20 x 16 in.	Desk	29 in.	1 ea/cs
CW0006CS	500 lb.	22 x 18 in.	Desk	30.5 in.	1 ea/cs
CW0006CSF	500 lb.	22 x 18 in.	Full length	30.5 in.	1 ea/cs
CW0007CS	500 lb.	24 x 18 in.	Desk	33 in.	1 ea/cs
CW0008CS	700 lb.	26 x 20 in.	Desk	34.5 in.	1 ea/cs
CW0009CS	700 lb.	30 x 20 in.	Desk	39 in.	1 ea/cs
Swing away footrests — powder coat finish					
CW0070PS	300 lb.	16 x 16 in.	Desk	25 in.	1 ea/cs
CW0071PS	300 lb.	18 x 16 in.	Desk	27 in.	1 ea/cs
CW0005PS	300 lb.	20 x 16 in.	Desk	29 in.	1 ea/cs
CW0006PS	500 lb.	22 x 18 in.	Desk	30.5 in.	1 ea/cs
CW0007PS	500 lb.	24 x 18 in.	Desk	33 in.	1 ea/cs
CW0008PS	700 lb.	26 x 20 in.	Desk	34.5 in.	1 ea/cs
CW0009PS	700 lb.	30 x 20 in.	Desk	39 in.	1 ea/cs
Elevating leg rests — chrome finish					
CW0070CEL	300 lb.	16 x 16 in.	Desk	25 in.	1 ea/cs
CW0071CEL	300 lb.	18 x 16 in.	Desk	27 in.	1 ea/cs
CW0071CELF	300 lb.	18 x 16 in.	Full length	27 in.	1 ea/cs
CW0005CEL	300 lb.	20 x 16 in.	Desk	29 in.	1 ea/cs
CW0006CEL	500 lb.	22 x 18 in.	Desk	30.5 in.	1 ea/cs
CW0006CELF	500 lb.	22 x 18 in.	Full length	30.5 in.	1 ea/cs
CW0007CEL	500 lb.	24 x 18 in.	Desk	33 in.	1 ea/cs
CW0008CEL	700 lb.	26 x 20 in.	Desk	34.5 in.	1 ea/cs
CW0009CEL	700 lb.	30 x 20 in.	Desk	39 in.	1 ea/cs
Elevating leg rests — powder coat finish					
CW0070PEL	300 lb.	16 x 16 in.	Desk	25 in.	1 ea/cs
CW0071PEL	300 lb.	18 x 16 in.	Desk	27 in.	1 ea/cs
CW0005PEL	300 lb.	20 x 16 in.	Desk	29 in.	1 ea/cs
CW0006PEL	500 lb.	22 x 18 in.	Desk	30.5 in.	1 ea/cs
CW0007PEL	500 lb.	24 x 18 in.	Desk	33 in.	1 ea/cs
CW0008PEL	700 lb.	26 x 20 in.	Desk	34.5 in.	1 ea/cs
CW0009PEL	700 lb.	30 x 20 in.	Desk	39 in.	1 ea/cs

With **over 30 different combinations available**, we have the wheelchair that is right for you.

Removable desk-length armrests

Wide, deep chair with high weight capacity

Elevating leg rests that swing away for easy patient transfer

Add an anti-tip device



CW0007CEL



CW0005CS



CW0005PEL



CW0006PS



CW00071PS



CW0005PS



CW0006CS



CW0007PEL



CWT0004CEL

Transport chairs

Cat. no.	Weight capacity	Dimensions (width x depth)	Description	Packaging
CWT0004CEL	300 lb.	19 x 16 in.	Chrome, elevating leg rest, full length armrest	1 ea/cs
CWT0004PEL	300 lb.	19 x 16 in.	Powder-coat, elevating leg rest, full length armrest	1 ea/cs

Wheelchairs and accessories (continued)

Wheelchair accessories

Cat. no.	Description	Packaging
CAX4028IV	Oxygen tank/IV pole	1 ea/cs
CAX20072	Oxygen tank holder	1 ea/bx; 3 bx/cs
CAIV20080	IV holder	1 ea/bx; 3 bx/cs
CAT20073	Lock assembly (for elevated or swing away foot rests)	1 pr/bg; 50 bg/cs
CATD20074	Anti-tip, rear	1 pr/bg; 10 bg/cs
CAAF20075	Anti-fold (for 300 lb. wheelchairs)	1 ea/bx; 20 bx/cs
CAAF20075B	Anti-fold (for 500 lb. wheelchairs)	1 ea/bx; 20 bx/cs
CAANT20014	Anti-theft device (fits 16 in. to 24 in. wheelchairs)	1 ea/bx; 2 bx/cs



CATD20074



CAIV20080

Wheelchair replacement parts

Cat. no.	Description	Packaging
CAAR20076	Desk length armrest	1 pr/bx; 5 pr/cs
CAAR20087	Full length armrest	1 pr/bx; 5 pr/cs
CAEL20077	Elevated leg rest	1 pr/bx; 5 pr/cs
CAS20078	Swing away foot rest	1 pr/bx; 5 pr/cs
CACP20019	Calf pad, black	1 pr/bg; 10 bg/cs
CAAP200DL	Desk arm pad, black	1 pr/bg; 10 bg/cs
CACST001	Front castor, 300 lb.	1 pr/bx; 5 pr/cs
CACST002	Front castor, 500 lb.	1 pr/bx; 5 pr/cs
CAUP200796	Replacement upholstery, 16 in. (width)	1 ea/bg; 20 bg/cs
CAUP200798	Replacement upholstery, 18 in. (width)	1 ea/bg; 20 bg/cs
CAUP200790	Replacement upholstery, 20 in. (width)	1 ea/bg; 20 bg/cs
CAUP200792	Replacement upholstery, 22 in. (width)	1 ea/bg; 20 bg/cs
CAUP200794	Replacement upholstery, 24 in. (width)	1 ea/bg; 20 bg/cs



CAAF20075



CAAR20076



CAAR20087

Note: The wheelchair accessories listed above are not compatible with the 700 pound weight capacity wheelchairs. Please contact your local Cardinal Health sales representative for more information.



CAEL20077



CAS20078



CAAP200DL



CACST001



CACST002



CAUP200790

Rollators and walkers

Choose from a wide range of weight-bearing options (for rollators: 250, 300 and 400 pounds—and for walkers: 300 and 500 pounds). Once again, you'll find what you need for bariatric patients right here.

Rollators

Cat. no.	Weight capacity	Dimensions (Handle height)	Dimensions (Seat height)	Description	Packaging
CRL0007	300 lb.	29–34 in.	20 in.	Rollator, black, junior, aluminum, 5 in. wheels	1 ea/cs
CRL0008	250 lb.	32–38 in.	20 in.	Rollator, blue, aluminum, 5 in. wheels	1 ea/cs
CRL0009B	400 lb.	35–40 in.	23.5 in	Rollator, black, heavy duty, steel, 5 in. wheels	1 ea/cs

Walkers

Cat. no.	Weight capacity	Dimensions (Walker height)	Width (Back legs)	Recommended user height	Description	Packaging
CWAL0006	300 lb.	25.5–32.5 in.	22 in.	4 ft. 8 in.–5 ft. 6 in. (Junior)	Walker, dual release, aluminum, single front 5 in. wheels	2 ea/cs
CWAL0008T	300 lb.	32.5–39.5 in.	23 in.	5 ft. 6 in.–6 ft. 4 in. (Adult)	Walker, dual release, aluminum, single front 5 in. wheels	2 ea/cs
CWAL0010B	500 lb.	33–43 in.	27 in.	5 ft. 6 in.–6 ft. 4 in. (Heavy duty adult)	Walker, dual release, aluminum, dual front 5 in. wheels	2 ea/cs
CWAL0007N	300 lb.	32–39 in.	23 in.	5 ft. 6 in.–6 ft. 4 in. (Adult)	Walker, dual release, aluminum, no wheels	2 ea/cs
CWAL1028H	500 lb.	29.5–37 in.	23 in.	5 ft. 1 in.–6 ft. (Adult)	Side hemi walker	2 ea/cs
CWAL240KS	300 lb.	Adjustable	-	One size fits most	Knee scooter, steerable, 8 in. wheels	1 ea/cs

Walker replacement parts and accessories

Cat. no.	Description	Packaging
CAWW1001	Single 5 in. wheels, pair	1 pr/bx; 25 bx/cs
CAWW2002	Dual 5 in. wheels, heavy duty, pair	1 pr/bx; 10 bx/cs
CAWW3003	Walker leg replacement, 8 holes	1 pr/bx, 25 bx/cs



CRL0009B



CRL0007



CWAL240KS



CWAL1028H



CWAL0006



CWAL0007N



CWAL0010B



CAWW1001

Crutches and replacement parts

After a few days on crutches, you know how uncomfortable they can get. Our answer for sore armpits?

Replacing the traditional foam rubber with **Thermoplastic Rubber (TPR)**—a synthetic blend that lasts long, resists peeling and most important, is soft. Your patients will feel the difference—and so will you, as you help them heal.

Crutches

Cat. no.	Weight capacity	Recommended user height	Description	Packaging
CA901AD	300 lb.	5 ft. 2 in.–5 ft. 10 in. (Adult)	Axillary, push button, adjustable, aluminum	8 pr/cs
CA801ADB	650 lb.	5 ft. 2 in.–5 ft. 10 in. (Adult)	Axillary, heavy duty, push button, adjustable, steel	4 pr/cs
CA901TL	300 lb.	5 ft. 10 in.–6 ft. 6 in. (Tall)	Axillary, push button, adjustable, aluminum	8 pr/cs
CA801TLB	650 lb.	5 ft. 10 in.–6 ft. 6 in. (Tall)	Axillary, heavy duty, push button, adjustable, steel	4 pr/cs
CA901YTH	300 lb.	4 ft. 6 in.–5 ft. 2 in. (Youth)	Axillary, push button, adjustable, aluminum	8 pr/cs
CA901CH	300 lb.	4 ft. 0 in.–4 ft. 6 in. (Child)	Axillary, push button, adjustable, aluminum	8 pr/cs
CA801FYTH	250 lb.	5 ft. 3 in.–6 ft. 3 in. (Tall)	Forearm, push button, adjustable, aluminum	10 pr/cs
CA801FCH	250 lb.	5 ft. 0 in.–6 ft. 1 in. (Adult)	Forearm, push button, adjustable, aluminum	10 pr/cs

Crutch replacement parts

Cat. no.	Description	Packaging
CACP20028	Underarm pad, gray	100 ea/cs
CACP20028B	Underarm pad, heavy duty, gray	100 ea/cs
CAACK20026	Crutch tip, underarm pad and hand grip	50 pr/cs
CAACK20026B	Parts kit, heavy duty, crutch tip, underarm pad and hand grip	50 pr/cs
CAACT20027	Crutch tip, gray	100 ea/cs
CAACT20027B	Crutch tip, heavy duty, gray	100 ea/cs



CA901AD



CA801ADB



CA801FCH




CAACK20026B

Bathroom accessories

Bathroom visits are more routine when you use the right products. Ours are recommended by clinicians, tested by experts and provide you with the highest quality at an affordable price. Cleaning instructions are simple and clear. We take care of the products, so you can do the same for your patients.

Dimension guide:



Commodes and accessories

Cat. no.	Weight capacity	Dimensions (Height)	Dimensions (Width)	Dimensions (Depth)	Description	Packaging
CPAT0002	300 lb.	17–21 in.	24 in.	16 in.	Commode, fixed arm, with lid, aluminum, includes 10 qt. bucket	4 ea/cs
CPAT0009	300 lb.	19.5–23.5 in.	24 in.	20 in.	Commode, drop arm, with lid, steel, powder-coat, includes 10 qt. bucket	2 ea/cs
CPAT0009D	300 lb.	19–23.5 in.	24 in.	20 in.	Commode, drop arm, with lid, aluminum, includes 7 qt. bucket	2 ea/cs
CPAT003F	300 lb.	16–22 in.	22 in.	20 in.	Commode, folding, with lid, 3-in-1, powder-coat, includes 7 qt. bucket	4 ea/cs
CPAT0008D	800 lb.	21 in.	26 in.	21 in.	Commode, bariatric, drop arm, steel, powder-coat, includes 7 qt. bucket	1 ea/cs
CABU0010	-	6.1 in.	11 in.	7 qt.	Bucket, 7 qt., with splash guard, with lid	10 ea/cs
CABU0010B	-	8.5 in.	11.25 in.	10 qt.	Bucket, heavy duty, 10 qt., with splash guard, with lid	10 ea/cs
CABU0007B	-	6 in.	11 in.	7 qt.	Bucket, square, 7 qt., with lid for use with CPAT0008D	10 ea/cs



CPAT0002



CPAT0009



CPAT0008D



CABU0010B

Toilet seats and accessories

Cat. no.	Weight capacity	Dimensions (Height)	Description	Packaging
CBAS0026	300 lb.	4.5 in.	Toilet seat, without clamp (fits most toilet seats)	1 ea/bx; 3 bx/cs
CBAS0027	300 lb.	5.5 in.	Toilet seat, with clamp (fits most toilet seats)	1 ea/bx; 3 bx/cs
CBAS0028	600 lb.	4.5 in.	Toilet seat, heavy duty, with legs and handles	1 ea/bx; 3 bx/cs
CSR00041	300 lb.	Adjustable	Toilet safety rail	1 ea/bx; 6 bx/cs



CBAS0026



CBAS0027



CBAS0028

Shower chairs, transfer bench and accessories

Cat. no.	Weight capacity	Dimensions (Seat height)	Description	Packaging
CBAS0029	250 lb.	14.3–20.3 in.	Shower chair, without back	2 ea/cs
CBAS0029T	250 lb.	14.3–20.3 in.	Shower chair, with back	2 ea/cs
CBAS0031	550 lb.	16.5–20.5 in.	Shower chair, heavy duty, with back	2 ea/cs
CBAS0032	375 lb.	21.25 in.	Commode shower chair, with back, with locking casters, includes 10 qt. bucket	2 ea/cs
CBAS0030	250 lb.	18.5–22.5 in.	Transfer bench, with back, aluminum	2 ea/cs
CACT2689	-	-	Shower chair tips, gray, rubber	4 ea/bg; 25 bg/cs



CBAS0029



CBAS0031



CBAS0032



CBAS0030



CACT2689

Canes, replacement parts and daily living aids

Canes

Cat. no.	Weight capacity	Dimensions (Height)	Base type	Handle type	Description	Packaging
CNE00013V	250 lb.	30–39 in.	Single	J-hook	Cane, push button, vinyl	25 ea/cs
CNE0012	250 lb.	30–39 in.	Single	J-hook	Cane, push button, aluminum, black	25 ea/cs
CNE0014	250 lb.	30–39 in.	Single	Offset	Cane, push button, aluminum, black	25 ea/cs
CNE0001	300 lb.	36.5 in.	Single	J-hook	Cane, natural wood	25 ea/cs
CNE0020B	500 lb.	19–38 in.	Single	Offset	Cane, heavy duty, push button, steel, black	6 ea/cs
CNE0021	300 lb.	30–39 in.	Quad	Offset	Cane, push button, aluminum, black	4 ea/cs
CNE0022B	500 lb.	29–38 in.	Quad	Offset	Cane, heavy duty, push button, steel, black	2 ea/cs



Cane replacement parts

Cat. no.	Description	Packaging
CAN20023Q	Cane tip, quad, black	100 ea/cs
CAN20023QB	Cane tip, heavy duty, quad, black	100 ea/cs
CAN20024	Cane tip, standard, black (17mm or 3/4 in.)	100 ea/cs
CAN20024W	Cane tip, wood cane, black	100 ea/cs



Daily living aids

Cat. no.	Description	Packaging
CRTX2009DE	Cart oxygen tank, D and E cylinder wheels	2 ea/cs



Products built with
patients in mind:

Wheelchairs and accessories

Transport chairs

Walkers and rollators

Crutches and replacement parts

Bathroom accessories

Canes, replacement parts and
daily living aids

Moving care forward

One source for low cost and high quality products

For more information or to order, contact your Cardinal Health sales representative today, call **877.CARDINAL** or visit **cardinalhealth.com/dme**.



cardinalhealth.com/dme