


**Your mask is an
important decision.
Choose the right one.**





Cardinal Health has simplified and standardized its facial protection product line to make your choice easier. We offer three lines of high-quality masks—each designed to address your requirements and budget.





Think about all the hazards in today's OR—
from blood to body fluids. Is the mask you're
wearing protecting you from these hazards?
Your mask is an important decision.
Choose the right one.



ASTM Level 3 Masks offer clinicians the protection they seek, without compromising the comfort and breathability they expect.

The four layers of ASTM Level 3 Masks

ASTM Level 3 Masks feature a proprietary four-layer construction, which offers you protection, breathability and comfort in the OR.



Cardinal Health™ ASTM Level 3 Masks meet the highest performance requirements defined by ASTM F2100-11

Medical face mask material requirements by performance level*


Characteristic	Level 3	Level 1
Bacterial filtration efficiency	≥98%	≥95%
Sub-micron particulates filtration efficient at 0.1 micron	≥98%	≥95%
Differential pressure, mm H2O/cm2 (Breathability)	<5.0	<4.0
Resistance to penetration by synthetic blood, minimum pressure in mm Hg for pass result	160 mm Hg	80 mm Hg
Flame spread	Class 1	Class 1

*Tested in accordance to ASTM F2100-11. Data on file.


Cardinal Health™ ASTM Level 3 Masks

ASTM Level 3 Surgical Masks

Not made with natural rubber latex



AT74535



AT74635

Cat. no.	Description	Qty. bx	Qty. cs
AT74535	Level 3 Surgical Mask with anti-fog foam	50	300
AT76835	Level 3 Surgical Mask with anti-fog adhesive	50	300
AT744235	Level 3 Surgical Mask with lightweight anti-fog strip	50	300
AT744335	Level 3 Surgical Mask with lightweight anti-fog strip and eyeshield	25	100
AT74635	Level 3 Surgical Mask with anti-fog foam and eyeshield	25	100
AT74635-I	Level 3 Surgical Mask with anti-fog foam and anti-glare eyeshield	25	100
AT54535	Level 3 Duckbill-Style Surgical Mask with anti-fog foam	50	300
AT54635	Level 3 Duckbill-Style Surgical Mask with anti-fog foam and eyeshield	25	100
AT54635-I	Level 3 Duckbill-Style Surgical Mask with anti-fog foam and anti-glare eyeshield	25	100

ASTM Level 3 Procedure Mask

Not made with natural rubber latex



AT74531



AT74631

Cat. no.	Description	Qty. bx	Qty. cs
AT74531	Level 3 Procedure Mask with anti-fog foam	50	300
AT74631	Level 3 Procedure Mask with anti-fog foam and eyeshield	25	100

Cardinal Health™ ASTM Level 1 Masks

ASTM Level 1 Masks are the standard for both surgical and procedural masks. As part of the complete Cardinal Health facial protection offering, they are a cost-effective solution for general use, particularly in low-fluid cases.

ASTM Level 1 Surgical Masks

Not made with natural rubber latex



Cat. no.	Description	Qty. bx	Qty. cs
AT71035	Level 1 Surgical Mask	50	300
AT71039	Level 1 Surgical Mask with horizontal ties	50	600
AT7505-P	Level 1 Surgical Mask (Poly inner layer)	50	300
AT71235	Level 1 Surgical Mask with anti-fog foam	50	300
AT51334	Level 1 Duckbill Surgical Mask (Poly inner layer)	50	300
AT51035	Level 1 Duckbill Surgical Mask (Cellulose inner layer)	50	300
AT73035	Level 1 Surgical Mask Dye-Free, white (Cellulose inner layer)	50	300
AT73335	Level 1 Surgical Mask Dye-Free, white (Poly inner layer)	50	300
AT752005	Level 1 Surgical Mask with anti-fog film	50	300
AT752007	Level 1 Surgical Mask with anti-fog film, XL	50	300
AT72835	Level 1 Surgical Mask with anti-fog adhesive	50	300
AT73835	Level 1 Surgical Mask with anti-fog odorless adhesive	50	300
AT75835	Level 1 Surgical Mask with anti-fog adhesive and removeable liner	50	300
AT771145A	Level 1 Surgical Mask with pediatric print	50	300


Cardinal Health™ ASTM Level 1 Masks

ASTM Level 1 Procedure Masks

Not made with natural rubber latex



AT7511



AT7004

Cat. no.	Description	Qty. bx	Qty. cs
AT7511	Level 1 Procedure Mask, blue (Poly inner layer)	50	300
AT7511-WE	Level 1 Procedure Mask with eyeshield, blue (Poly inner layer)	25	100
AT71021	Level 1 Procedure Mask, blue (Tissue inner layer)	50	500
AT70021	Level 1 Procedure Mask, yellow (Tissue inner layer)	50	500
AT7004	Level 1 Procedure Mask, yellow (Poly inner layer)	50	300
AT771141	Level 1 Adult Size Earloop Procedure Mask with pediatric print	50	300

Specialty Masks

Not made with natural rubber latex



AT771212



AT7509†

Cat. no.	Description	Qty. bx	Qty. cs
AT771212	Pediatric Face Cover	50	750
AT7509†	Cone Mask	50	300

†Cone mask is not a high-filtration mask.

Cardinal Health™ Surgical N95 Respirators

Our cost-effective Surgical N95 Respirator line meets Centers for Disease Control and Prevention (CDC) guidelines for protection against TB and is National Institute for Occupational Safety and Health (NIOSH) certified to have a filter efficiency level of 95 percent or greater against particulate aerosols free of oil. These Surgical N95 Respirators are also Food and Drug Administration (FDA) cleared for use in the OR during surgical procedures.

N95 respirators

Not made with natural rubber latex



USA-N95-S



N95-ML

Cat. no.	Description	Qty. bx	Qty. cs
USA-N95-S	Flat-fold Surgical N95 Respirator, small	50	300
USA-N95-R	Flat-fold Surgical N95 Respirator, regular	50	300
N95A-S	Molded Surgical N95 Respirator, small	20	120
N95A-ML	Molded Surgical N95 Respirator, medium/large	20	120

Cardinal Health™ Visor Systems

Cardinal Health offers a comprehensive portfolio of protective eyewear to help meet your eyewear needs.

The Cardinal Health™ Protective Eyewear System includes:



Eyewear

- A flexible, lightweight frame
- Optical-grade, clear replacement lenses
- Pre-assembled eyewear frame and lens option*

Cat. no.	Qty. bx	Qty. cs
CAHFRAME	n/a	100
CAHLENS	100	600
CAHEYEWEAR [†]	n/a	50

Shields

- One-piece face shield
- Foam headband
- $\frac{3}{4}$, full-length and extended-length options

Cat. no.	Qty. bx	Qty. cs
H1SHIELD50 ($\frac{3}{4}$ length)	25	50
F1SHIELD50 (full-length)	25	50
XLSHIELD50 (extended-length)	20	50

[†]Note: CAHEYEWEAR consists of a pre-assembled CAHFRAME and CAHLENS



Comprehensive solutions to support the
entire continuum of care

For those tasked with navigating the complexities of healthcare, Cardinal Health brings scaled solutions that help our customers thrive in a changing world. We offer a broad portfolio of products and services that span across the continuum of care — from hospital to home.

We are Essential to Care™



To help educate your staff on choosing the right mask for the right procedure, please access our free facial protection CE offering: cardinalhealth.com/education



Facial Protection