



Cardinal Health™ Monoject™ Magellan™ Insulin and Tuberculin Safety Syringes

Clinically preferred

Safety needles designed to protect



Your facility depends on solutions that offer greater protection against needlesticks and bloodborne pathogens

Because more than 40% of needlestick injuries happen after use,¹ we designed Cardinal Health™ Monoject™ Magellan™ Safety Needles to keep patients and clinicians safe at every step. Magellan™ Safety Needles are engineered to comply with OSHA regulations and help you standardize solutions across every job with one safety technology.

Choose Monoject™ Magellan™ Insulin and Tuberculin Safety Syringes

Magellan™ Products are the clinically preferred premium option for safety needles and syringes.

Intuitive, easy-to-use designs offer protection with just the push of a finger, improving clinician comfort and confidence with less follow-up training.

- 1 Robust safety shield with low profile design provides coverage of entire cannula
- 2 One-handed activation, two ways
- 3 Keeps hand behind needle
- 4 Audible activation cue
- 5 Minimal in-servicing required

Every Monoject™ Magellan™ Insulin and Tuberculin Safety Syringe does not require a change in protocol and may be activated using the thumb or finger.



Cardinal Health™ Monoject™ Magellan™ Insulin and Tuberculin Safety Syringes

Features:

- Sterile, single use and individually packaged
- Supports diabetes care
- Not manufactured with DEHP or natural rubber latex



Monoject™ Magellan™ Insulin Safety Syringes with Permanent Needles

Order code	Description	Syringe volume	Needle gauge	Needle length	Box qty.	Case qty.
8881892910	1 mL 29 G x ½ in. (0.33 x 13 mm)	1 mL	29 G	½ in. (13 mm)	50	500
8881892950	½ mL 29 G x ½ in. (0.33 x 13 mm)	½ mL	29 G	½ in. (13 mm)	50	500
8881892930	¾ mL 29 G x ½ in. (0.33 x 13 mm)	¾ mL	29 G	½ in. (13 mm)	50	500
8881893010	1 mL 30 G x ⅝ in. (0.3 x 8 mm)	1 mL	30 G	⅝ in. (8 mm)	50	500
8881893050	½ mL 30 G x ⅝ in. (0.3 x 8 mm)	½ mL	30 G	⅝ in. (8 mm)	50	500
8881893030	¾ mL 30 G x ⅝ in. (0.3 x 8 mm)	¾ mL	30 G	⅝ in. (8 mm)	50	500

Note: ¾ mL = 30 units, ½ mL = 50 units, 1 mL = 100 units

Monoject™ Magellan™ Tuberculin Safety Syringes with Permanent Needles

Order code	Description	Syringe volume	Needle gauge	Needle length	Box qty.	Case qty.
8881882712	1 mL 27 G x ½ in. (0.4 x 13 mm)	1 mL	27 G	½ in. (13 mm)	50	500
8881882812	1 mL 28 G x ½ in. (0.36 x 13 mm)	1 mL	28 G	½ in. (13 mm)	50	500

Monoject™ Magellan™ Tuberculin Safety Syringe Trays

Order code	Description	Syringe volume	Needle gauge	Needle length	Box qty.	Case qty.
8881882712T	1 mL 27 G x ½ in. (0.4 x 13 mm)	1 mL	27 G	½ in. (13 mm)	20 each/tray	40 trays/case
8881882812T	1 mL 28 G x ½ in. (0.36 x 13 mm)	1 mL	28 G	½ in. (13 mm)	20 each/tray	40 trays/case



For more information or to order Monoject™ Sharps Safety Products, contact your sales representative, call **800.964.5227** or visit **cardinalhealth.com/monoject**

Access the Monoject™ Magellan™ Insulin and Tuberculin Safety Syringes instructions for use at **cardinalhealth.com/monoject/IFU**

References: 1. NIOSH. Alert: preventing needlestick injuries in health care settings. Cincinnati, OH: US Dept of Health and Human Services; 1999. <https://www.cdc.gov/niosh/docs/2000-108/pdfs/2000-108.pdf>.

© 2021 Cardinal Health. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO, MONOJECT and MAGELLAN are trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other trademarks are the property of their respective owners. Patent cardinalhealth.com/patents. Lit. No. 2GM21-1652238-01 (09/2021)

