



Cardinal Health™ Monoject™ Veterinary Needle and Syringe Selection Guide



Monoject™ Veterinary Aluminum Hub Needles and Syringes




Monoject™ Veterinary Metal Insert Needles and Syringes



Monoject™ Hospital Standard Needles and Syringes



Monoject™ Magellan™ Safety Needles and Syringes



The choice is yours with Cardinal Health™ Monoject™ Needles and Syringes

Choose your needle

- Veterinary angle point: Designed for strength to penetrate thick or dehydrated skin
- Hospital point: Designed for smooth penetration into thin, well-hydrated skin
- Metal insert for robust needle attachment
- Epoxy needle attachment
- Clear or colored hub and sheath

Choose your syringe

- Luer lock or slip tip syringe
- Unique silver or white plunger

Choose your packaging

- Recyclable rigid pack
- Peel-to-open soft pack*

Cardinal Heath™ Monoject™ Veterinary Needle and Syringe Selection Guide

Product features		Monoject™ Veterinary Aluminum Hub Needles and Syringes	Monoject™ Veterinary Metal Insert Needles and Syringes	Monoject™ Hospital Standard Needles and Syringes	Monoject™ Magellan™ Safety Needles and Syringes
Needle sheath	Color coding	•	•		
	Clear			•	•
Needle point	Shallow needle angle point for strength to penetrate thick or dehydrated skin	•	•		
	Steep needle angle point for penetration into thin, hydrated skin			•	•
Needle hub	Aluminum*	•			
	Clear stabilizing metal insert		•		
	Color-coded epoxy			•	•
Syringe attachment	Slip tip	•	•		
	Luer lock	•	•	•	•
Syringe graduations	Heat etched	•	•	•	
	Printed	•		•	•
Finger flanges	Low profile	•	•	•	
	Wide finger	•		•	•
Plunger color	Silver	•	•	•	
	White	•		•	•
Packaging	Recyclable rigid packaging	•	•	•	
	Peel-to-open soft packaging	•		•	•

“A survey of 200 animal health professionals found that over 92% of vet clinics reported having more than one needlestick injury per year.”*

– Blinded Veterinarian Study, January 2021 (data on file)



Consider using a safety engineered alternative device.

Contact your distributor sales representative to order.

For more information call **877.CARDINAL** or visit **cardinalhealth.com/veterinarysupplies**.



Reference:

*Data on file. Blinded Veterinarian Study, January 2021.