

See the Difference. Feel the Quality.

Uni-Patch™ C Series (First Choice) Stimulating Electrodes


Our Uni-Patch™ C series (first choice) stimulating electrodes are designed for durability, flexibility, and comfort, with:

- A long-lasting 0.04-inch thick hydrogel for reliable adhesion that conforms well over body contours
- A quick, moderately aggressive tack for firm yet gentle electrode-to-skin contact for multiple applications
- Highly conductive carbon film for even current distribution
- Spun lace top coat or soft foam covers for flexible and comfortable application to various sites
- Flexible, durable, and highly conductive pigtail wire strands designed for multiple uses

The soft mold flex connector with secure lead wire adhesive technology features:

- A socket that easily adapts to variations in pin size and provides a secure connection
- Wire adhesive designed to anchor the lead wire into the pad and diminish separation

When you're done using the electrodes, you can return them to either side of the double-coated, tinted storage liner.



Cloth stimulating electrodes

ORDERING INFORMATION

Uni-Patch™ C Series (First Choice) Stimulating Electrodes

Cloth stimulating electrodes

CODE	MODEL	SIZE	SHAPE	QTY
EP85795	3100C	1.25 in.	Round	4/pack
EP85765	3105C	2 in.	Round	4/pack
EP85770	3110C	2.75 in.	Round	4/pack
EP85775	3115C	2 in.	Square	4/pack
EP85780	3120C	2 x 3.5 in.	Rectangle	4/pack

10–14 applications

Uni-Patch™ C Series (First Choice) Stimulating Electrodes

Foam stimulating electrodes

CODE	MODEL	SIZE	SHAPE	QTY
EP85800	3150F	1.25 in.	Round	4/pack
EP85805	3155F	2 in.	Round	4/pack
EP85825	3160F	2.75 in.	Round	4/pack
EP85835	3165F	2 in.	Square	4/pack
EP85840	3170F	2 x 3.5 in.	Rectangle	4/pack

10–14 applications

For more information, please contact your Cardinal Health sales representative at OEMSpecialty@cardinalhealth.com or visit cardinalhealth.com/pain-therapy-recovery.

© 2018 Cardinal Health. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO and Uni-Patch are trademarks or registered trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other marks are the property of their respective owners. 02/2018 - H9829 - WF#755787

15 Hampshire Street
Mansfield, MA 02048

800.962.9888
508.261.8000

