

OrthoWedge™ Healing Shoe

Effective Off-Loading Footwear




The DARCO OrthoWedge™ is clinically proven to reduce weight bearing pressure on the forefoot which promotes faster healing after surgery, trauma or when forefoot wounds or ulcerations are present.

DARCO




OrthoWedge™ Features and Benefits:

- › **15° Wedge Sole** effectively shifts body weight to the midfoot and heel. Reduces forefoot pressure by **as much as 57%!**
- › **Square Toe Design** acts as a bumper to protect the toes when k-wires are present and provides better Left/Right fit.
- › **Ankle Strap** seats foot firmly in shoe to reduce heel slippage and eliminate friction.
- › **Removable Forefoot Closure** provides security while eliminating buckle pressure.
- › **Zoned Outsole** has an aggressive tread under the mid-foot where traction is needed most.
- › **Removable Insole** features twice the padding of standard insoles and can be modified as needed. Insole can be totally removed to accommodate DARCO's innovative **Peg Assist™** customizable off-loading insole.



Standard Surgical Shoe



DARCO OrthoWedge™

OrthoWedge™ Indications:

- › Post Surgical healing
- › Wounds or ulcerations present under metatarsal heads and toes.
- › Rearfoot trauma
- › Any situation that requires a pressure off-load of the forefoot.

Item Number	Description	Size	Color	Women's Size	Men's Size
OQ0B	OrthoWedge	X-Small	Black	4 - 7	
OQ1B	OrthoWedge	Small	Black	7.5 - 10	6 - 8
OQ2B	OrthoWedge	Medium	Black	10.5 - 13	8.5 - 10
OQ3B	OrthoWedge	Large	Black		10.5 - 12
OQ4B	OrthoWedge	X-Large	Black		12.5 - 14

Available through:



Headquarters:

DARCO International | 810 Memorial Blvd. | Huntington, WV 25701 | USA
 Toll Free voice 800-999-8866 | Fax 877-201-7082 | Voice 304-522-4883 | Fax 304-522-0037
 info@darcointernational.com | www.darcointernational.com

As with all DARCO products, a portion of the profits from the sale of these items will be donated to encourage and support continuing education in our profession.