


Axess™ Low Profile Nasal Mask

- Low-profile design
- Fits snugly around nose without tape
- Translucent
- Not made with natural rubber latex
- Available in Birthday Bubblegum scent, Fresh Mint scent, and unscented



Also available: Axess™ Capnography Adapter; attaches between the Axess™ Nasal Mask and the Axess™ Scavenging Circuit

Access to Hard-to-Reach Areas

The low-profile design of the Axess™ Nasal Mask provides unencumbered mouth access while its small, lightweight scavenging system is easy to work around.

Access to Greater Patient Comfort

Available in three sizes and scents, Axess™ Nasal Masks fit snugly without tape or uncomfortable tubes protruding into the patient's nostril to maximize patient comfort.

Access to Excellent Care – at a Lower Cost

Each Axess™ Nasal Mask comes individually wrapped for greater infection control, and comes at a cost that's lower than that of other competing single-use disposable nitrous oxide/oxygen masks.




"It really does smell like bubblegum!"

Nitrous Oxide/Oxygen Sedation Solutions

Comfort, Relaxation, Efficiency

The reusable gas delivery and scavenging system ensures that the desired amount of gas is delivered and expelled efficiently for your practice—and comfortably for your patients.



Axess™ Scavenging Circuits

Ref. #	Description
52006	Scavenging Circuit for RFS™ Remote Flow System and Accu-Vac™ Controllers
52005	Scavenging Circuit For In-Line Controller

Special
Introductory
Kits

Axess™ Low Profile Nasal Mask

Ref. #	Description	Quantity
53037-9	Small Clear Unscented	24/Box
53037-12	Small Birthday Bubblegum	24/Box
53037-16	Small Fresh Mint	24/Box
53035-9	Medium Clear Unscented	24/Box
53035-12	Medium Birthday Bubblegum	24/Box
53035-16	Medium Fresh Mint	24/Box
53034-9	Large Clear Unscented	24/Box
53034-12	Large Birthday Bubblegum	24/Box
53034-16	Large Fresh Mint	24/Box

Axess™ Capnography Adapters

Ref. #	Description	Quantity
53055	Single-use Capnography Adapters	24/Box
53004	Scavenging Circuit Extension	

Order from your preferred dealer today.

Learn more:



[Crosstex.com/Axess](https://www.crosstex.com/Axess)



All product names are trademarks of Crosstex International, Inc., a Cantel Medical Company, its affiliates or related companies, unless otherwise noted. Marks not registered in all jurisdictions.

© 2020 Crosstex International, Inc. DLIT000019 Rev C 0120