





HCPCS:

- Control panel activates and adjusts water pressure and cleansing modes
- Bidet provides front or female and back cleansing depending on needs
- Retractable bidet nozzle and guard for a more sanitary experience
- Contoured for safe and comfortable elevation
- Antimicrobial* treated raised toilet seat helps minimize the growth of bacterial odors, mold, and fungi on treated surfaces.
- *Antimicrobial properties are built-in to protect the products. These products do not protect users or others against bacteria, viruses, germs, or other disease organisms.
- RTL12C005-WH Ambient water temperature only
- RTL12C006-WH Choose between warm and ambient water temperatures
- RTL12C005-WH Easily connects to the water source behind the toilet
- RTL12C006-WH Easily connects to the water source behind the toilet and under the sink
- RTL12C006-WH Comes with an inlet cap to plug warm water connection for ambient water use only to accommodate bathroom layouts

Description
RTS w. Bidet, Ambient Water, 1/CV
RTS w. Bidet, Warm Water, 1/CV
Description
One Year Limited
HDPE
400 lbs
White
E0244
no