

## Pure Expressions Dual Channel Electric Breast Pumps

HCPCS: E0603-1,E0603



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- LCD Display provides status of programs and functions
- Independent & custom suction for each breast 9 adjustable suction strengths
- 2 Phase Expression Stimulation and Expression Mode
- Includes 2 Soft Silicone Cushion Inserts for Optimum Comfort
- Easy to assemble and clean Silent for discreet pumping
- Closed Loop Anti-backflow System
- Single or Double Pumping Options
- BPA and Latex Free
- Two 4oz Bottle Feeding Adapter Kits and Storage Containers with Lids
- Convenient Carry Bag with Zipper
- Insulated Cooler Bag comes with RTLBP2200 only
- Includes A/C Power Supply & External battery pack 4 AA alkaline batteries (not included)
- One Year Limited Warranty

PRODUCTS	Description
RTLBP2000	Double Electric Plus Breast Pump 1/ea
SPECIFICATIONS	Description
Batteries (pair)	4 AA Alkaline
Contains Latex	No
Product Dimensions	4" (W) x 6.25" (L) x 2.75" (H)
Product Weight	.9 lbs
Retail Box Dimensions	RGLBP200 12.2" (W) x 7.5" (L) x 11" (H) RTLBP2200 12.2" (W) x 9.1" (L) x 12.8" (H)
Retail Box Weight	5.4 lbs.
Warranty	One Year Limited
ACCESSORIES	Description
BP0015	FeedingNippleSet, BP1000/2000/2200,1/pk, PKG
BP0016	Base for Bottles, BP1000/2000/2200, 1/e, EA
BP0013	PumpKit w/o Insert,31MM,BP1000/2000/220, BAG
BP0011	PowerAdapter,6V,RTLBP2000/RTLBP2200, 1e, EA
BP0012	ExternalBatteryPack,6V,RTLBP2000/BP2200, EA
BP0014	Silicone Insert, 27MM, BP1000/2000/2200, BAG
BP008	Membranes ,BP1000/2000/2200, 4/bag, BAG
BP009	Tubing Set, BP1000/2000/2200, 1/pr, PR
BP006	Silicone Insert, 24MM, BP1000/2000/2200, BAG

## ACCESSORIES

## Description

BP007	ValveSetW/Membrane,BP1000/2000/2200,1/p, PR
BP004	6oz wide neck bottle w/cap. 1/ea, EA
BP005	PumpKit w/o Insert,27MM,BP1000/2000/220, BAG
BP002	Insulated Cooler Bag, 1/ea, EA
BP001	Breast Pump Carry Bag,1/ea, EA